WINCARE

**Axtair  
Automorpho<sup>®</sup>  
Plus**

Monographie

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<sup>1</sup> **LOCF**: Last Observation Carried Forward

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10 FÉVRIER 2021

[winnicare.com](http://winnicare.com)

# Incidence of pressure ulcer in patients using an alternating pressure mattress overlay: the ACTIVE study

**Objective:** The primary objective was to determine the clinical benefit of using a specific alternating-pressure mattress overlay (APMO) in the prevention of pressure ulcer (PU) in patients at medium to high risk.

**Method:** This prospective study was conducted in five rehabilitation centres and three nursing homes. Patients at medium to high risk of PU, but without PU at baseline, and lying between 15 and 20 hours per day on a specific APMO were included. The primary endpoint was the percentage of patients who developed a sacral, spine, heel or trochanteric PU (supine support areas) of at least category II, at day 35. All patients were included in the analysis.

**Results:** A total of 89 patients were included; of whom six patients (6.7%) dropped out of the study (average ( $\pm$ standard deviation) follow-up 32 $\pm$ 5.4 days). No sacral, spine, heel or trochanteric PU of at least category II was reported (i.e., an incidence of 0% [95%

Confidence Interval: 0–4.1%] according to the exact Clopper–Pearson method]. Patients were ‘satisfied’ or ‘very satisfied’ with the comfort and stability of the APMO. The caregivers assessed as ‘very easy’ or ‘easy’ the implementation, maintenance and use of the APMO (turning over, moving to a sitting position).

**Conclusion:** In combination with the usual measures to prevent PU, the results of our study showed a low incidence of PU in high-risk patients lying for between 15 and 20 hours a day on an APMO, use of which is therefore recommended in these patients.

**Declaration of interest:** SM received honoraria as coordinator of the study indirectly by Nukleus from Winnicare France, France. CR and MM (working for Nukleus) also received indirect compensation from Winnicare France. The study was sponsored and funded by Winnicare France but it did not participate in the conduct of the study nor in the analysis of the data.

alternating-pressure mattress overlay • beds • life support system • pressure ulcer • prevention • wound • wound care • wound healing

A pressure ulcer (PU) (also called pressure injury, pressure sore or bedsore) is a localised injury to the skin and/or underlying tissue, usually over a bony prominence, as a result of pressure, or pressure in combination with shear.<sup>1,2</sup>

The lesion may be in the form of intact skin or an open wound and may be painful. Ischaemia, stress, recovery of blood flow, tissue hypoxia and the pathological impact of pressure and shear pathological mechanisms interact in the development of PU. Soft tissue tolerance to pressure and shear forces can also be affected by microclimate, nutrition, tissue perfusion, comorbidities, and the condition of the skin and underlying tissue.<sup>1</sup>

The prevalence and burden of PU remain high,<sup>3</sup> and understanding of PU pathophysiology is progressing.<sup>4,5</sup>

PU prevention requires a global approach. It has been the subject of recommendations.<sup>1,6–8</sup> Among the measures to prevent PU in patients who are bedbound for a long time every day, it is recommended to reduce the duration of pressure and its intensity on the bone areas.<sup>1</sup> Several types of support surface designed to reduce the intensity and duration of pressure on contact points are available to caregivers in order to prevent PU. In 2018, the National Pressure Ulcer Advisory Panel Support

Surface Standards Initiative<sup>9</sup> specified the foundational definition of a support surface: ‘A specialized device for pressure redistribution designed for management of tissue loads, microclimate, and/or other therapeutic functions (e.g., any mattresses, integrated bed system, mattress replacement, overlay, or seat cushion, or seat cushion overlay).’ The NPUAP has classified support surfaces as basic/standard hospital mattress and reactive support surface (non-powered or powered).<sup>9</sup>

The technical specifications of each support surface determine their expected performance in terms of PU prevention. For alternating-pressure air mattresses or mattress overlays, air cell thickness, flow and inflation/deflation cycle, and physical properties play a very important role in their effectiveness in terms of prevention.<sup>10</sup>

The present study aimed to determine the clinical benefit of using a specific alternating-pressure mattress overlay (APMO) in the prevention of pressure ulcer (PU)

Sylvie Meaume,<sup>1</sup> MD, Head of Service\*; Chantal Rousseaux,<sup>2</sup> MD, Responsible for Medical Affairs; Marc Marty,<sup>1</sup> MD, Methodologist

\*Corresponding author email: m.marty@nukleus.fr

<sup>1</sup> Rothschild Hospital – Geriatric Department and Wound Care Unit, Assistance Publique Hôpitaux de Paris, Paris, France. <sup>2</sup> Nukleus, Clinical Research Department, Paris, France.

in patients at medium to high risk. The specific characteristics of the mattress are described in the method section.

The risk of PU was reduced by 69% in the alternating-pressure dynamic air mattress/mattress overlay groups versus standard mattress in the 2015 meta-analysis of McInnes et al.<sup>10</sup> (risk ratio (RR)=0.31; 95% confidence interval (CI): 0.17–0.58) including two trials. The risk of methodological bias was considered high. The difference between alternating-pressure devices versus constant low-pressure devices was not established.<sup>10</sup>

In a network meta-analysis, powered active air surfaces showed a statistical superiority in PU prevention in comparison with standard mattresses with an estimated RR of 0.42 (CI: 0.29–0.63) (moderate certainty of evidence) and with non-powered reactive foam surfaces with an estimated RR of 0.64 (CI: 0.42–0.96) (low certainty of evidence).<sup>11</sup> According to a recent systematic review, it is necessary to continue investigations to assess the use of support surfaces in the prevention of PU.<sup>12</sup>

The choice of surface support must be made according to each patient's level of risk of PU occurrence.<sup>13,14</sup> The use of each support surface must be adapted to the level of risk of PU occurrence and to the context.

The conditions for reimbursement by French national health insurance for support surfaces to assist in the prevention and treatment of PUs were specified in the 'Commission nationale d'évaluation des dispositifs médicaux et des technologies de santé' (CNEDIMTS) opinion of December 2009.<sup>15</sup> The CNEDIMTS is an independent commission issuing opinions for health regulatory authorities in France. For alternating-pressure surface supports, manufacturers must provide evidence attesting to technical specifications and clinical data attesting to the performance of the support surface. Clinical data must follow the recommendations of the CNEDIMTS, appearing in Annex III to the opinion of December 2009.<sup>15</sup> According to the specifications of the mattresses or mattress overlay, indications for a prevention and adjuvant device in the treatment of PUs are precisely specified.

The main objective of this study was to determine the clinical benefit of using a specific APMO in the prevention of PU in patients at medium to high risk of PU.

## Method

### Design

This was a national, prospective, multicentre, observational, longitudinal study, with a patient monitoring/follow-up period of 35 days ( $\pm 5$  days). The study took place between June and December 2019 in five rehabilitation centres and in three nursing homes. Investigators were physicians specialised in physical medicine and rehabilitation, and neurologists or geriatricians qualified to take care of patients at risk of PU. The study was non-interventional and did not involve any risk or constraint in any of the procedures performed, and the products were used in the usual way, according to clinical guidelines.

The protocol follows the recommendations of the CNEDIMTS appearing in Annex III to the opinion of December 2009.<sup>15</sup>

### Patients

Patients had to meet the following criteria to participate in the study: >18 years old; at medium to high risk of PU (clinical judgement and a score between 10–14 on the Braden scale<sup>16</sup> (six, maximum risk to 23, no risk); without PU on the day of inclusion; woken up during the day; lying between 15–20 hours per day on a specific APMO (but for <48 hours on this APMO); and weighing between 30–165kg. Patients (or their representative) should have been informed of the study and agreed to participate. In the event of the patient's inability to give consent, the patient's person of trust/legal representative could give consent for the patient to take part in the study. To fulfil CNEDIMTS requirements, patients with a life expectancy of <6 months or with malnutrition could not be included.

Malnutrition was defined for adults <70 years of age as weight loss  $\geq 5\%$  in one month or  $\geq 10\%$  in six months or a body mass index (BMI)  $\leq 18.5\text{kg/m}^2$  (excluding constitutional weakness); for adults  $\geq 70$  years of age, weight loss  $\geq 5\%$  in one month or  $\geq 10\%$  in six months or BMI  $\leq 21\text{kg/m}^2$  or the Mini Nutritional Assessment (MNA)<sup>17</sup> score  $\leq 17$  (/30) or albumin  $< 35\text{g/l}$ .<sup>18</sup>

### Course of study and data collection

When a patient fulfilled the participation criteria, the evaluating physician offered the possibility of participating in the study to the patient or to their representative, in the event of the patient's incapacity. If consent was given to participate, an information leaflet was given to the patient or their representative. A baseline visit was made, and the patient then continued to be monitored by the care team, which used the standard PU prevention measures, such as measures to decrease pressure (changing position, installation, use of supports), observing skin condition, maintaining skin hygiene, ensuring nutritional balance, promoting patient/caregiver participation in PU prevention and ensuring continuity of care.

Skin condition was monitored daily to detect the slightest appearance of a PU throughout the 35-day follow-up period. On day 35 ( $\pm 5$  days), a final visit was planned to record the occurrence of PUs during the study and their progression in the case of occurrence, and the secondary criteria of evaluation.

At baseline, the following information was recorded: patient demographics (age, sex, weight, height); place of patient care; condition responsible for the risk of PU; Braden score; comorbidities; previous mattress support surface (before installation of the specific APMO used in the study); skin condition—absence of current PU and history of PU; average time spent lying per day; bed installation; seated installation (equipment/support); physiotherapy; activity level of the patient; nursing protocol implemented (number of position changes per day; maximum time between two position changes; number of sheet changes per day, positions used).

At the end of the study (day 35 or before in the case of withdrawal), the occurrence of one or more PUs was recorded and, where applicable, their characteristics (location, size, worst category during the period, progression). The category was considered according to the NPUAP classification.<sup>19</sup> The following parameters were also recorded: patient's view on the comfort of the mattress (general comfort and stability); satisfaction of the nursing staff (ease of implementation, maintenance, use in terms of turnaround, in terms of moving the patient to a sitting position) (on a scale of 0–4); degree of moisture (humidity component of the Braden score); sound level of the mattress (on a scale of 0–4). If a patient had difficulty in responding on the comfort of the mattress because of their neurological state, the experience of the patient was recorded by the care staff in consultation with the patient's representative.

Technical incidents related to the mattress (such as breakdowns, valve problems) and adverse events (AE) were also recorded. An information log book (paper) was provided to patients and caregivers to collect any important clinical events.

The protocol specified that, in the event of premature discontinuation of the study of a patient, the patient should be evaluated in a similar way as on day 35.

#### **Pressure-relieving support and measures to prevent PU**

To be included in this observational study, patients had to lie for <48 hours on the APMO (Axtair Automorpho Plus, Winnicare France, France). This APMO comprises an overlay, a foam support and a compressor, and is available in three bed widths (90cm, 100cm and 120cm).

The overlay comprises 18 independent and removable polyurethane ether cells with a 12cm height of therapeutic air, a head area with two static cells covered with a pillow, and four heel cells with independent discharge. The overlay is laid on a support base in polyether foam, >5cm thick, located in an independent compartment. A cardiopulmonary resuscitation (CPR) valve was added to permit a mattress deflation time of <20 seconds, assuming the need for emergency cardiopulmonary resuscitation. The device (overlay and base) is fully protected by a removable cover made of material impermeable to liquids but permeable to water vapour.

The main mode of action of the mattress is active. The air cells of the upper layer (overlay) are inflated and deflated regularly by the compressor. The compressor is based on a patented system of automatic and continuous calculation of the inflation pressure according to the patient's size. There are two other programmable operating modes available: static low-pressure mode, allowing the treatment of patients requiring transient immobilisation, and care mode, allowing the mattress overlay to be overinflated for up to 30 minutes to facilitate transfers or treatments. In practice, the use of these modes over the period of a day are short and intended for the actions of caregivers. Visual and audible alarms were available, and the operating manual was provided to patients and caregivers.

A previous study was carried out on 57 patients at medium or high risk of PU (Braden score  $\leq 17$ , bedbound for >15 hours per day), free from PU at baseline in rehabilitation centres and using this specific APMO. Over a follow-up period of 31 days, a PU incidence of 7.0% (95% CI: 2.0–17.0) was highlighted.<sup>20</sup> All patients received PU prevention measures as personalised care protocols for each patient.

#### **Outcomes**

The primary endpoint was the percentage of patients who developed a sacral, spine, heel or trochanteric PU (supine support areas) between day 0 and day 35 or the end of the study. PU was defined as skin damage of at least stage/category 2 of the NPUAP classification.<sup>19</sup> Daily skin inspection made it possible to detect the occurrence of any PU over the follow-up period. Stage/category I PUs (persistent redness of the skin >24 hours) were not taken into account in the primary endpoint because they are difficult to diagnose with certainty (differential diagnosis problems) and because of large inter-individual variation in their assessment.<sup>10</sup>

The secondary end points were category I (all areas) PU incidence or PU incidence of any category in an area other than sacrum, heel, spine, heels and trochanteric; patient satisfaction with the comfort of the APMO; patient acceptance of the sound level of the APMO; and the care team's assessment of the use of the APMO and the moisture level. Safety was analysed by description of AEs and technical incidents.

#### **Sample size**

Assuming that 7% of patients developed a PU of at least category II of the sacrum, spine, heel and trochanteric between day 0 and day 35, 80 patients were needed to make it possible to demonstrate that the upper limit of the 95% CI of the percentage would be <20%, with a power of >95%. In view of the data in the literature,<sup>21–26</sup> a maximum upper bound of 20% appears to be clinically relevant.

#### **Analysis of data**

All the data from patients who agreed to participate in the study were taken into account, according to the general principle of intention-to-treat analysis.

No statistical test was carried out, the study being of a descriptive nature. Any missing data were extrapolated using the last observation carried forward technique. This meant that, in the event of premature termination of the study, the patient had to be evaluated and the data for the last day the patient was in the study was used as the data for day 35.

The 95% CI for the percentage of patients who developed a PU between inclusion and day 35 was calculated using the exact Clopper–Pearson method. Secondary endpoints are reported descriptively.

#### **Ethical and regulatory approval**

The project was submitted to an Ethics Committee (Personal Protection Committee Northwest I - Rouen

France No. 2019-A00798-49) and was approved May 2019.

Protection and confidentiality of patient data were guaranteed by compliance with applicable laws.

## Results

### Enrolment and baseline characteristics

A total of 89 patients were included in the study, of whom six (6.7%) dropped out of the study. Reasons include:

- Patient could not bear the inflation of the tubes (n=1)
- Patients presented a non-device-related AE (n=2, one patient with pulmonary congestion requiring hospitalisation on day 15, and one patient with an abscess on a scar following surgery on the cervical spine and requiring hospitalisation on day 14)
- Patient needed a change of mattress on day 17 (problem of bed width for obesity) (n=1)
- Patient needed a change of mattress on day 21 (unknown reason) (n=1)
- Patient returned home (n=1).

The average time of follow-up (time between day zero and the last evaluation) was  $32 \pm 5.4$  (standard deviation, SD) days for the 89 patients. Patients (of whom 57.3% were male) were on average  $73.1 \pm 20.5$  years old and had an average BMI of  $25.3 \pm 6.2 \text{ kg/m}^2$ . Demographic and clinical characteristics of patients at baseline are shown in Table 1. Of the patients, 38 (42.7%) stayed in five rehabilitation centres and 51 (57.3%) patients in three nursing homes. The average length of stay for patients was  $2.6 \pm 3.7$  years (n=89). Urinary incontinence was present in 85.4% of patients and anal incontinence in 68.5% of patients.

Patients included were at high risk of PU—mean score on the Braden scale was  $12.8 \pm 1.6$  at inclusion—with an average moisture score (on a scale of 0–4) of  $2.2 \pm 0.8$ . Patients were lying down on average  $16.6 \pm 1.8$  hours/day. No patient had a PU at inclusion and 18 (20.2%) had a history of PU. The medical conditions responsible for the risk of PU were classified into: neurological (for example, hemiplegia, tetraplegia, stroke) for 44 (49.4%) patients; accidental (for example, fracture, trauma) for 25 (28.1%) patients; or multifactorial or of complex origin (such as cancer, older age) for 16 (17.9%) patients. In bed, 52 (58.4%) patients used cushions to prevent PUs, 18 (20.2%) patients used positioning materials and four (4.5%) patients used orthoses. When seated, 69 patients (77.5%) used a cushion to prevent PU.

### Primary end point

No sacral, spine, heel or trochanteric PU of at least category II was reported (i.e., an incidence of 0%; 95% CI: 0–4.1%, according to the exact Clopper–Pearson method). The upper limit of the 95% CI was well below the 20% assumption and attests to the performance of the APMO.

### Secondary end points: performance

A category I sacral PU occurred nine days after installation on the APMO in an 85-year-old female

patient with dementia, with urinary and anal incontinence, a Braden score of 10 at baseline, and who was bedbound for 15 hours/day. The PU was still progressing at the end of the study on day 32 (no further data were available after this day) (without worsening). No PU in any other location was reported.

Patients were ‘satisfied’ or ‘very satisfied’ with the APMO in the majority of cases in terms of comfort (87.5%) and stability (77.3%) (Table 2).

The care team assessed as ‘very easy’ or ‘easy’ in the majority of cases the implementation (93.2%) and maintenance of the APMO (95.5%), and its use in terms of turning over (85.2%) and sitting position (88.6%) (Table 3).

Patients rated the sound level of the APMO as satisfactory or very satisfactory in the majority of cases (93.1%) (Table 2).

At day 35, the degree of moisture had not changed as compared to day 0. The moisture score was 2.2 (SD 0.8) on day 0 and 2.3 (SD 0.9) on day 35 on the moisture component of the Braden scale (ranging from 1 (constantly moist) to 4 (rarely moist)), and did not change between baseline and the end of the study.

### Secondary endpoints: adverse events

In total, three AEs were reported in three patients, including two patients who stopped the study prematurely:

- An 85-year-old patient who developed a category I sacral PU nine days after installation on the study mattress
- A 46-year-old male patient, weighing 104kg and measuring 176cm, with quadriplegia of level C3 with left hemiplegia, bedbound 15 hours a day, prematurely released from the study for hospitalisation on day 15 due to pulmonary congestion. This event was considered a serious AE unrelated to the APMO. The patient did not develop PUs
- A 53-year-old male patient with accident-related mobility problems, bedbound 20 hours a day, weighing 60kg and measuring 170cm, hospitalised on day 14 for an abscess on a scar. This event was considered a serious AE unrelated to the APMO. The patient did not develop PUs.

### Secondary endpoints: technical incidents

A technical incident was reported: in one case there was a difficulty in sliding the patient onto the shower stretcher, due to the mattress cover hanging over the mattress.

## Discussion

The number of PUs observed during this study is very low, with an incidence of 0% (95% CI: 0–4.1% according to the exact Clopper–Pearson method) PUs in the included population. The upper limit of the 95% CI was well below the 20% assumption.

These results on PU incidence are similar to those observed in a similar study conducted in patients at risk of PUs and lying 10–15 hours per day on a similar

APMO.<sup>26</sup> The PU incidence observed in the MATCARP study was 1.2% (1/83).<sup>26</sup> Indirect comparisons should always be interpreted with caution because of the many confounding factors that can affect the results, such as the residence of the patients, their degree of disability, the course of disability (acute or chronic stage), PU prevention methods and evaluation methods.<sup>27</sup> Nevertheless, these results should be considered, bearing in mind the recent results of the very large PRESSURE<sup>2</sup> randomised clinical trial conducted in 2029 patients in 42 UK secondary/community facilities. The trial compared alternating-pressure mattresses (APMs) with high specification foam (HSF) for a maximum treatment phase of 60 days.<sup>28</sup> The primary outcome was time to developing a new PU category  $\geq$ II from randomisation to 30 days from the end of the treatment phase (maximum 90 days). The incidence of at least one new category II PU (7.9%) was not statistically significantly different between the two groups, with an absolute difference of 2% (APMs 70 PUs (6.9%), HSF 90 PUs (8.9%)). In a treatment phase sensitivity analysis (until day 90), 132 (6.5%) patients developed a new category II PU between randomisation and the end of treatment phase (APM 53 (5.2%), HSF 79 (7.8%)), with a statistically significant difference observed in time to development of category II PU in a Fine and Gray model (hazard ratio, HR=0.66 (95% CI: 0.46–0.93; exact p=0.0176)).<sup>28</sup>

The choice of patients included in the study was based on the criteria retained by the French authorities for this type of APMO (in this case, an alternating-pressure mattress or APMO comprising between 10–15cm of therapeutic air) to allow access to reimbursement, provided that the support device guarantees precise minimum technical specifications. The French authorities justify the exclusion of malnourished patients and patients with a life expectancy of <6 months to avoid the issue of people being lost to follow-up. The choice of patient profile at risk of PU by the French authorities in 2009 and corresponding to who should be treated with this type of APMO was largely based on expert advice. Indeed, the choice of a surface support adapted to the risk of the patient to prevent PUs remains a question that cannot be completely based on a high level of evidence.<sup>13</sup>

### Limitations

This observational study presents certain weaknesses which attenuate the generalisation of the results. The non-comparative nature of the study alters the causality between the clinical outcome and the use of the APMO. Indeed, other preventive PU measures were used for each patient. These consisted of limiting the pressure time (change of position, installation), regularly observing skin condition, maintaining skin hygiene, ensuring nutritional balance, promoting the participation of the patient and their caregiver, and ensuring the continuity of care. The causality of the low incidence of PU observed in the study is therefore

**Table 1. Baseline characteristics of patients (full analysis set, n=89). Figures are mean±standard deviation (SD) unless specified**

	Patients, n=89
Age, year	73.1±20.5
Male gender, n (%)	51 (57.3)
Body mass index, kg/m <sup>2</sup>	25.3±6.2
Urinary incontinence, n (%)	76 (85.4)
<b>Urinary incontinence (n=76), n (%)*</b>	
Intermittent	23 (31.1)
Total	51 (68.9)
Anal incontinence, n (%)	61 (68.5)
<b>Anal incontinence (n=61), n (%)</b>	
Intermittent	15 (24.6)
Total	46 (75.4)
Braden score (6–23)	12.8±1.6
Hours in bed per day	16.6±1.8
<b>Installation in bed, n (%)</b>	
Use of cushions	52 (58.4)
Use of orthotics	4 (4.5)
Medical device for patient positioning	18 (20.2)
Number of changes of position per day (n=31)	5.2±1.4
Time between two changes of position (n=30), hours	4.0±1.2
Number of linen changes per day (n=32)	1.0±0.3
Using pressure-relieving cushions when seated, n (%)	69 (77.5)
<b>Condition responsible for the risk of pressure ulcers, n (%)</b>	
Injury	25 (28.1)
Multifactorial	16 (17.9)
Neurological	44 (49.4)
Other	20 (22.5)
History of pressure ulcers, n (%)	18 (20.2)
*Data on the nature of the urinary incontinence is missing for two patients	

limited. Nevertheless, in view of the high risk level of occurrence of PU in the patients included, the use of this type of support surface appears to be suitable and it seems unlikely that the APMO has not contributed to PU prevention.

The duration of the study was only 35 days and may be considered too short. Even if it is established that a PU can appear in a few minutes or hours, a longer study period than 35 days could have shown a higher incidence. This period of >30 days is sufficient for the French authorities for this type of study, which is why we have chosen it.

**Table 2. Patients' opinion on the comfort and sound level of the APMO at end of the study\***

	General comfort n=88 <sup>†</sup>	Stability n=88 <sup>†</sup>	Sound level n=87 <sup>‡</sup>
	n, (%)		
Not at all satisfactory	2 (2.3)	2 (2.3)	—
Not satisfactory	3 (3.4)	8 (9.1)	1 (1.2)
Neither satisfactory nor unsatisfactory	6 (6.8)	10 (11.4)	5 (5.8)
Satisfactory	51 (58.0)	44 (50.0)	21 (24.1)
Very satisfactory	26 (30.0)	24 (27.3)	60 (69.0)

APMO—alternating-pressure mattress overlay; <sup>†</sup>For 17 patients, due to neurological state, patient experience was recorded by the staff with subjective assessment by caregivers and/or patient representative(s); <sup>‡</sup>one missing data; <sup>§</sup>two missing data

**Table 3. Care teams' opinions on using the APMO, n=88\***

	Set-up	Ease of cleaning	Ease of use in terms of patient turnaround	Ease of use in terms of lying to sitting
	n, (%)			
Not at all satisfactory	—	—	—	—
Not satisfactory	—	—	—	—
Neither satisfactory nor unsatisfactory	6 (6.8)	4 (4.6)	13 (14.8)	10 (11.4)
Satisfactory	26 (29.5)	25 (28.4)	55 (62.5)	31 (35.2)
Very satisfactory	56 (63.6)	59 (67.1)	20 (22.7)	47 (53.4)

APMO—alternating-pressure mattress overlay; \*one missing data

Malnourished and end-of-life patients were excluded from the study. This requirement of the French authorities for this type of study limits the generalisation of the results, malnourished and end-of-life patients being particularly exposed to the risk of PU. This requirement was made to avoid patients leaving studies prematurely.

However, the study has strengths. All of the patients were taken into account in the analysis and there was no attrition bias. The patients studied are well described and appear to be appropriate for the type of APMO used. STROBE recommendations were followed to optimise the quality of data presentation.<sup>29</sup>

These clinical data, which complement technical specifications about the support surface, provide some evidence of the performance of the APMO studied in prevention of PU.

Patients were satisfied or very satisfied with the comfort and stability of the APMO in the majority of cases. The care teams assessed as very easy or easy in the majority

of cases the implementation and maintenance of the APMO, and its use in terms of turning over and moving to a sitting position. The patients rated the sound level of the APMO as 'satisfactory' or 'very satisfactory' in the majority of cases. For 17 patients, the patient experiences were recorded by the staff and the family of the patient by subjective assessment, and so these criteria should be analysed with caution. No APMO-related AE was reported. Only one technical incident occurred (difficulty sliding the patient onto the shower stretcher, the mattress cover often hanging on the mattress). All of these data attest to the performance of the APMO studied.

**Conclusion**

In conclusion, in combination with the standard PU prevention measures, the study data indicated a low incidence of PU in the medium to high risk patients in this study, who were lying for between 15 and 20 hours a day on an APMO. **JWC**

**Reflective questions**

- On what criteria do you currently use an alternating pressure mattress overlay (APMO) for patients at risk of ressure ulcer (PU)?
- How might the results of this study change your opinion on the use of an APMO in patients at risk for PU?
- How might the results of this study change your practice on the use of an APMO in patients at risk for PU?

## References

- 1 National Pressure Ulcer Advisory Panel, European Pressure Ulcer Advisory Panel and Pan Pacific Pressure Injury Alliance. Prevention and treatment of pressure ulcers: quick reference guide. Emily Haesler (Ed.). Cambridge Media, 2014
- 2 National Pressure Ulcer Advisory Panel (NPUAP) announces a change in terminology from pressure ulcer to pressure injury and updates the stages of pressure injury. 2016. [tinyurl.com/vpt4ozp0](https://tinyurl.com/vpt4ozp0) (accessed 28 January 2021)
- 3 Moore Z, Avsar P, Conaty L et al. The prevalence of pressure ulcers in Europe, what does the European data tell us: a systematic review. *J Wound Care* 2019; 28(11):710–719. <https://doi.org/10.12968/jowc.2019.28.11.710>
- 4 Mervis JS, Phillips TJ. Pressure ulcers: pathophysiology, epidemiology, risk factors, and presentation. *J Am Acad Dermatol* 2019; 81(4):881–890. <https://doi.org/10.1016/j.jaad.2018.12.069>
- 5 Blackburn J, Ousey K, Taylor L et al. The relationship between common risk factors and the pathology of pressure ulcer development: a systematic review. *J Wound Care* 2020; 29(Sup3):S4–S12. <https://doi.org/10.12968/jowc.2020.29.sup3.s4>
- 6 National Institute for Health and Care Excellence. Pressure ulcers: prevention and management. Clinical guidelines CG179. <https://www.nice.org.uk/guidance/cg179> (accessed 28 January 2021)
- 7 Agency of Health Care Research and Quality. Pressure ulcer risk assessment and prevention: comparative effectiveness. 2013. [tinyurl.com/y8ouc4u6](https://tinyurl.com/y8ouc4u6) (accessed 28 January 2021)
- 8 Qaseem A, Mir TP, Starkey M, Denberg TD. Risk assessment and prevention of pressure ulcers: a clinical practice guideline from the American College of Physicians. *Ann Intern Med* 2015; 162(5):359. <https://doi.org/10.7326/M14-1567>
- 9 National Pressure Ulcer Advisory Panel. Support surface standards initiative terms and definitions related to support surfaces. Version 01/29/2007 Update 12/27/2018. [tinyurl.com/yhzr6tgg](https://tinyurl.com/yhzr6tgg) (accessed 2 February 2021)
- 10 McInnes E, Jammali-Blasi A, Bell-Syer SE et al. Support surfaces for pressure ulcer prevention. *Cochrane Database Syst Rev* 2015; 9(9):CD001735. <https://doi.org/10.1002/14651858.CD001735.pub5>
- 11 Shi S, Dumville JC, Cullum N. Support surfaces for pressure ulcer prevention: a network meta-analysis. *PLoS ONE* 2018; 13(2): e0192707. <https://doi.org/10.1371/journal.pone.0192707>
- 12 Rae KE, Isabel S, Upton D. Support surfaces for the treatment and prevention of pressure ulcers: a systematic literature review. *J Wound Care* 2018; 27(8):467–474. <https://doi.org/10.12968/jowc.2018.27.8.467>
- 13 McNichol L, Mackey D, Watts C, Zuecca N. Choosing a support surface for pressure injury prevention and treatment. *Nursing* 2020; 50(2):1. <https://doi.org/10.1097/01.NURSE.0000651620.87023.d5>
- 14 McNichol L, Watts C, Mackey D et al. Identifying the right surface for the right patient at the right time: generation and content validation of an algorithm for support surface selection. *J Wound Ostomy Continence Nurs* 2015; 42(1):19–37. <https://doi.org/10.1097/WON.000000000000103>
- 15 Haute Autorité de Santé. Opinion of 22 December 2009 of the National Commission for the Evaluation of Medical Devices and Health Technologies on 'Medical devices for the prevention of pressure ulcers' (Title I, Chapter 2, section 1, sub-section 2 of List of products and services mentioned in Article L165-1 of the Social Security Code). [In French.] 2009. [tinyurl.com/265m84e2](https://tinyurl.com/265m84e2) (accessed 28 January 2021)
- 16 Braden BJ, Bergstrom N. Predictive validity of the Braden scale for pressure sore risk in a nursing home population. *Res Nurs Health* 1994; 17(6):459–470. <https://doi.org/10.1002/nur.4770170609>
- 17 Bauer JM, Kaiser MJ, Anthony P et al. The Mini Nutritional Assessment—its history, today's practice, and future perspectives. *Nutr Clin Pract* 2008; 23(4):388–396. <https://doi.org/10.1177/0884533608321132>
- 18 Opinion of the commission - Commission evaluation of products and benefits - September 27, 2006 Products: Home nutrition products and related services. <https://www.hassante.fr/portail/upload/docs/application/pdf/cepp-469.pdf> [AQ9: Cannot find reference, please provide new ref/link]
- 19 Edsberg LE, Black JM, Goldberg M et al. Revised national pressure ulcer advisory panel pressure injury staging system. *J Wound Ostomy Continence Nurs* 2016; 43(6):585–597. <https://doi.org/10.1097/WON.0000000000000281>
- 20 Barrois B, Gaubert-Dahan ML, Allaert FA et al. Prevention of pressure ulcers with a motorized air support in at-risk patients hospitalized in rehabilitation departments APAM + an observational prevalence study with historical control and incidence study. *Panminerva Med* 2018; 60(4):161–169. <https://doi.org/10.23736/S0031-0808.18.03491-2>
- 21 Nixon J, Cranny G, Iglesias C et al. Randomised, controlled trial of alternating pressure mattresses compared with alternating pressure overlays for the prevention of pressure ulcers: PRESSURE (pressure relieving support surfaces) trial. *BMJ* 2006; 332(7555):1413. <https://doi.org/10.1136/bmj.38849.478299.7C>
- 22 Vanderwee K, Grypdonck M, Defloor T. Alternating pressure air mattresses as prevention for pressure ulcers: a literature review. *Int J Nurs Stud* 2008; 45(5):784–801. <https://doi.org/10.1016/j.ijnurstu.2007.07.003>
- 23 Sanada H, Sugama J, Matsui Y et al. Randomised controlled trial to evaluate a new double-layer air-cell overlay for elderly patients requiring head elevation. *J Tissue Viability* 2003; 13(3):112–121. [https://doi.org/10.1016/S0965-206X\(03\)80016-5](https://doi.org/10.1016/S0965-206X(03)80016-5)
- 24 Serraes B, Beeckman D. Static air support surfaces to prevent pressure injuries: a multicenter cohort study in Belgian nursing homes. *J Wound Ostomy Continence Nurs* 2016; 43(4):375–378. <https://doi.org/10.1097/WON.0000000000000244>
- 25 Zwaenepoel E, Van Hecke A, Manderlier B et al. Pressure ulcer Cat. II-IV incidence on the CuroCell S.A.M. PRO powered reactive air support surface in a high-risk population: a multicentre cohort study in 12 Belgian nursing homes. *Int Wound J* 2020; 17(1):124–131. <https://doi.org/10.1111/iwj.13242>
- 26 Meaume S, Marty M. Pressure ulcer prevention using an alternating-pressure mattress overlay: the MATCARP project. *J Wound Care* 2018; 27(8):488–494. <https://doi.org/10.12968/jowc.2018.27.8.488>
- 27 Berlowitz D. Incidence and prevalence of pressure ulcers. In: Thomas DR, Compton G (eds). *Pressure ulcers in the aging population: a guide for clinicians*. Springer Science+Business Media New York; 2014: 19–26
- 28 Nixon J, Smith IL, Brown S et al. Pressure Relieving Support Surfaces for Pressure Ulcer Prevention (PRESSURE 2): clinical and health economic results of a randomised controlled trial. *EclinicalMedicine* 2019; 14:42–52. <https://doi.org/10.1016/j.eclinm.2019.07.018>
- 29 von Elm E, Altman DG, Egger M, et al.; STROBE Initiative. The Strengthening of Reporting of Observational Studies in Epidemiology (STROBE) statement: guidelines for reporting observational studies. *Ann Intern Med* 2007; 147(8):573–577. <https://doi.org/10.7326/0003-4819-147-8-200710160-00010>

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Conflicts of interest.—Brigitte Barrois, Marie-Line Gaubert-Dahan and François-André Allaert received grant to coordinate the study; Ilinca Urzica, Bernard Bordet, Jean-Marc Michel, Marie-Ange Blanchon and Cédric Chol received fees to be investigator of the study.

Acknowledgments.—The authors thank the director and paramedical staff of the participating nursing homes in Rehabilitation Services centers: Stéphanie Piot, Stéphanie Grocolas (CRR F COS Pasteur – Troyes); Dominique Delandre, Nikolay Dobrev, Mohammed Yahla, Laurence Leveque, Christelle Tardi, Céline Picard (CR F Pasori – Cosne-sur-Loire); Patricia Blaise, Adeline Rabczynski, Marc Allenbach, Claire Matter, Sophie Mavel (Les Hôpitaux Civils de Colmar – Colmar); Sylvie Favier, Jean-Christophe Poncet (Hôpital Charité du CHU de Saint-Etienne – Saint-Etienne). In addition, the authors thank Winnicare France (manufacturer) and these service providers (1er Mai médical, Bernamont médical, LVM88, HMS42) for providing the medical services needed for this study in accordance with good practice recommendations.

Manuscript accepted: June 4, 2018. - Manuscript received: May 31, 2018.

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[winnocare.com](http://winnocare.com)

ORIGINAL ARTICLE

# Prevention of pressure ulcers with a motorized air support in at-risk patients hospitalized in rehabilitation departments

## APAM + an observational prevalence study with historical control and incidence study

Brigitte BARROIS <sup>1</sup>, Marie-Line GAUBERT-DAHAN <sup>2</sup>, François-André ALLAERT <sup>3</sup>\*, Ilinca URZICA <sup>4</sup>, Bernard BORDET <sup>5</sup>, Jean-Marc MICHEL <sup>6</sup>, Marie-Ange BLANCHON <sup>7</sup>, Cédric CHOL <sup>7</sup>

<sup>1</sup>Department of Quality, Gonesse Hospital, Gonesse, France; <sup>2</sup>Department of Rehabilitation, Raincy-Montfermeil Hospital, Montfermeil, France; <sup>3</sup>Medical Evaluation Chair BSB and CEN Biotech, Dijon, France; <sup>4</sup>Department Rehabilitation, COS Pasteur, Troyes, France; <sup>5</sup>Department of Rehabilitation, Cosne-sur-Loire, France; <sup>6</sup>Department of Rehabilitation, Hospital of Colmar, Colmar, France; <sup>7</sup>Department of Rehabilitation, University Hospital, Saint-Etienne, France

\*Corresponding author: François-André Allaert, Medical Evaluation Chair BSB and CEN Biotech, Dijon, France. E-mail: [allaert@groupeccen.com](mailto:allaert@groupeccen.com)

### ABSTRACT

**BACKGROUND:** The main study was a cross-sectional multicenter study of the prevalence of pressure ulcers in rehabilitation services (RS) where Axtair Automorpho® Plus mattresses were made available *ad libitum*. The primary objective was to enable comparison with prevalence as observed in the national PERSE study (historical control) in similar departments using a variety of prevention aids.

**METHODS:** This cross-sectional prevalence study was combined with a prospective study of the incidence of pressure ulcer occurrence in the high-risk target population whose beds were systematically equipped with the studied support.

**RESULTS:** On the day of the prevalence study, 18 of the 456 patients presented pressure ulcers which occurred during the period when Axtair Automorpho® Plus supports were made available to the departments. This corresponds to a prevalence of 3.9% [2.4, 6.2] which is significantly lower than those recorded in the PERSE study: 11.8% [10.8; 12.8]  $p < 0.0001$ . The ulcer incidence study covers 57 patients who were hospitalized on a bed with an Axtair Automorpho® Plus support because of their risk of contracting pressure ulcers. One pressure ulcer occurred in 3 patients and 3 pressure ulcers occurred in another patient, i.e. 4 out of 57 corresponding to an incidence of 7.0% [2.0; 17.0].

**CONCLUSIONS:** The results of the pressure ulcer prevalence study in rehabilitation departments where Axtair Automorpho® Plus supports were available as well as results of the study of pressure ulcer incidence in patients bedridden on Axtair Automorpho® Plus, confirm the expected benefit provided to patients in terms of prevention.

(Cite this article as: Barrois B, Gaubert-Dahan ML, Allaert FA, Urzica I, Bordet B, Michel JM, *et al.* Prevention of pressure ulcers with a motorized air support in at-risk patients hospitalized in rehabilitation departments. Panminerva Med 2018;60:000-000. DOI: 10.23736/S0031-0808.18.03491-2)

**KEY WORDS:** Pressure ulcers - Incidence - Prevalence - Prevention- Axtair Automorpho® Plus.

Pressure ulcers are lesions of the skin or the subcutaneous tissues located over bony prominences resulting from the application of pressure, or pressure combined with shear forces.<sup>1</sup>

For more than 20 years, the PERSE Study has been conducting major national surveys in France of the prevalence of pressure ulcers and these have become benchmarks in the field.<sup>2-5</sup> The most recent survey conducted in 2015<sup>6</sup> up-

dated the data on the national prevalence of pressure ulcers, *i.e.* for all hospital departments and by type of department. It covered 776 hospital departments and 21,538 patients, 12,752 (59.2%) of whom were women and 8786 (40.8%) men. The survey shows that the prevalence of pressure ulcers was 8.1%, 95% CI = [7.7; 8.5] (N.=1753). This figure breaks down to 7.8%, 95% CI = [7.3; 8.3] (N.=997) in women and 8.6%, 95% CI = [8.0; 9.2] (N.=756) in men

( $P=0.0381$ ). In rehabilitation services, this prevalence was 11.8%, 95% CI = [10.8; 12.8]. Reducing the occurrence of pressure ulcers is a public health issue with respect to avoiding a major cause of deterioration in patients' quality of life and cutting the healthcare expenses they induce (2.7). In addition to the cost of care that they require in terms of medical supplies and nursing time, their occurrence is responsible for lengthening the average hospital stay by an estimated 9.8 days.<sup>7-10</sup> Pressure ulcers also require nursing care when patients return home, which in many places is not easy to arrange due to the shortage of nurses.

The onset of pressure ulcers can be prevented by identifying the risk factors presented by patients and by reducing the interface pressure between the skin and supporting surfaces.<sup>2</sup> Prolonged bed rest should be avoided by different methods, including patient mobilization, the use of different positions for bed and chair rest, the elimination of shear forces, and the use of items for preventing and treating pressure ulcers such as cushions, mattresses, and mattress toppers. As part of this prevention, the use of pressure relief mattresses is one of the five best practices with a major impact on risk management. Motorized air supports are considered to be of prophylactic and therapeutic value for patients at risk of developing moderate to severe pressure ulcers. Axtair Automorpho® Plus is a motorized air support consisting of a network of 18 air cells, 12 cm high, arranged on a 5 cm foam base and connected to a compressor. The compressor inflates the cells (one out of two) alternately and automatically regulates the inflation pressure. The alternation of pressures and the system for relieving pressure on the heels applied to four cells helps to prevent the prolonged vascular compression that may lead to tissue hypoxia causing the onset of pressure ulcers. To evaluate the contribution to reducing pressure ulcer prevalence and incidence from using this equipment in daily practice, we conducted an observational study with historical control in at-risk patients in rehabilitation services. The study method was approved beforehand by experts from France's national health evaluation authorities. This paper presents the results of the study which ran from 1 February to 30 June 2017.

## Materials and methods

### Experimental plan and primary objective

The main study was a cross-sectional multicenter study of the prevalence of pressure ulcers in rehabilitation ser-

vices (RS) where Axtair Automorpho® Plus mattresses were made available *ad libitum*. The primary objective was to enable comparison with prevalence as observed in the national PERSE Study (historical control) in similar departments using a variety of prevention aids. This cross-sectional prevalence study was combined with a prospective study of the incidence of pressure ulcer occurrence in the high-risk target population whose beds were systematically equipped with Axtair Automorpho® Plus mattresses.

### Secondary objectives

The secondary objectives were to describe the evolution of pre-existing pressure ulcers and of those occurring in patients during the study period, the satisfaction of the nursing staff, patient comfort, and the safe use of the support.

### Scientific committee of the study

Prior to the study, a scientific committee was appointed specifically to advise WINNCARE FRANCE on any medical or scientific issues arising during the study, on the exploitation of data, the drafting of the study report, and the publication of the results.

### Selection of RS centers

The study was conducted in rehabilitation services so that the population would be structurally similar to the elderly population maintained at home. The size of the sample required for the historical comparison meant that these centers had to be large enough to conduct a prevalence study on 330 patients and recruitment was conducted competitively in the various centers according to their size. All departments participating in the study were required to have existing pressure ulcer prevention nursing protocols which remained unchanged during the study.

### Patient inclusion criteria

In every rehabilitation service all incoming patients were to be included in the study except those who, after being informed, declined or could not consent because they were under 18 years of age or had severely impaired cognitive functions (MMSE Score <15), or because of a medical contraindication to the studied support. All patients were described in a screening log. The prevention equipment under study was Axtair Automorpho® Plus and patient beds were fitted with this equipment in line with the medical practitioners' decisions as and when patients were included in the study.

### Evaluation criteria

The primary endpoint was the occurrence or not of a pressure ulcer as defined by stages 1 to 4 of the NPUAP (National Pressure Ulcer Advisory Panel). This criterion was analyzed in terms of prevalence in the general population of hospitalized patients and in terms of incidence in the population of patients without pressure ulcers upon inclusion and meeting the severity criteria defined by the health authorities.

The secondary criteria were:

- the characteristics of pre-existing pressure ulcers and of ulcers that occurred in patients who were free of them when first hospitalized in the department;
- nursing staff satisfaction with the noise level of the prevention support, the ease of maintenance, and the ease of use in terms of turning patients and sitting them up;
- patient comfort and patients' opinions about the noise level of the equipment;
- product safety and in particular any adverse events and technical incidents with the product.

### Conduct of the study and data collection

#### *General organization of the study*

The study included a cross-sectional prevalence survey of all patients in services on the 60<sup>th</sup> day and 120<sup>th</sup> day after the start of the study plus a prospective survey of incidence over 5 months.

A register of all patients entering the services was maintained. It described their demographic characteristics, the reason for their stay, their risk factors for pressure ulcers, their skin condition and in particular the presence of pressure ulcers (location and description), patient's mobility, the assessment of the risk of pressure ulcers on the Braden scale, any contraindication to using Axtair Automorpho® Plus, whether or not the selection criteria were met, the patient's non-opposition to participating in the incidence study, and on completing the registration form, the decision whether or not to include the patient in the incidence study.

#### *Conduct of the cross-sectional study: prevalence assessment*

The prevalence study was conducted on all patients in the department on the 60<sup>th</sup> day and 120<sup>th</sup> day after the start of the study. The 60-day threshold was set because the average length of stay in these services was 31 days. By taking 60 days, save very few exceptions, all patients meeting the inclusion criteria had had the opportunity to benefit

from Axtair Automorpho® Plus prevention if the caregivers had deemed it necessary. For the same reason, it was very likely that by the 120<sup>th</sup> day all the patients in beds equipped with Axtair Automorpho® Plus would be different from those included in the 60-day prevalence analysis.

A specific observational register was completed for the assessments of the prevalence study. It described the number of men and women hospitalized, and the number of men and women with at least one pressure ulcer, regardless of its stage, and included a form for each patient describing their age, sex, reason for hospitalization, risk factors of pressure ulcers, pressure ulcer prevention supports, existence of pressure ulcers and their number, plus the stage and location of the most severe of them.

#### *Conduct of the prospective cohort study: incidence evaluation*

After verifying the patient's non-opposition to the collection of their personal information and validation of their selection criteria, the practitioner filled out an inclusion sheet on D0 describing the patient's demographic characteristics, the existence of any pressure ulcers and their description: location, length, width, severity, and clinical stage. The practitioner also took a standardized picture where applicable. For the evaluation of the incidence, the cutaneous state of the patients was reported on D15 and on D31, and also in the event of the occurrence of a pressure ulcer the day of its occurrence. Any pressure ulcers were described and photographed according to the same protocol as at the inclusion visit. The standardized photographic protocol required the picture to be taken 20 cm from the pressure ulcer with a colorimetric and centimeter scale. If the patient left the hospital before D31, a form was completed on the day of discharge.

A record was made at each visit, or in the event of transfer or discharge from the service before D31, of:

- an assessment of the risk of occurrence of the pressure ulcer (Braden scale);
- the maceration of the skin on a Likert scale from 0 to 5;
- unexpected events and possible technical incidents related to the Axtair Automorpho® Plus support;
- the evolution of any pressure ulcer that appeared;
- the patient's opinion about the comfort and noise level of the Axtair Automorpho® Plus support.

#### *Evaluation of the opinion of caregivers and patients*

On D15, the investigator recorded on a special form the caregivers' opinion of the Axtair Automorpho® Plus mo-

torized air support and their level of satisfaction with its use: noise level, ease of implementation and maintenance, ease of use in terms of turning patients and getting them to sit up evaluated on 5-point Lickert scales. Patient opinion and satisfaction with the comfort and noise level of Axtair Automorpho® Plus were also reported.

#### *Evaluation of tolerance*

Any adverse events were recorded on the follow-up questionnaire completed by the practitioners.

#### **Statistical analysis**

##### *Populations studied*

In each rehabilitation service, all patients were included in the study except those who, after being informed, declined or could not consent because they were under 18 years of age, or had severely impaired cognitive functions (MMSE<15), or a medical contraindication.

Populations for the prevalence study included:

- all the patients present in the services at the end of the 60<sup>th</sup> and 120<sup>th</sup> day, and presenting a moderate to high risk of pressure ulcers defined by a Braden score of 17 or less,<sup>11, 12</sup> who were up out of bed during the day, bedridden more than 15 hours, and who did not have any pressure ulcers on admission to the service;
- and all patients hospitalized in services present at the end of the 60<sup>th</sup> and 120<sup>th</sup> day.

Populations for the incidence study included:

- all patients hospitalized in departments whose beds were equipped with an Axtair Automorpho® Plus, having a hospital stay of at least 29 days, presenting a moderate to high risk of pressure ulcers defined by a Braden score of 17 or less, bedridden more than 15 hours, and not presenting any pressure ulcers on admission to the service;
- and all patients hospitalized in the service meeting the inclusion criteria whose beds were equipped with an Axtair Automorpho® Plus.

#### *Description of the statistical analysis*

##### DESCRIPTION OF THE POPULATIONS

Demographic and clinical characteristics of patients included in the prevalence and incidence populations were described by means and standard deviations for the quantitative variables and by numbers and percentages for the qualitative variables.

##### PREVALENCE AND INCIDENCE OF PRESSURE ULCERS

The prevalence of pressure ulcers in the rehabilitation services was described by number, frequencies, and 95% confidence intervals and compared to the current reference rate of the 2015 PERSE Study<sup>6</sup> by a  $\chi^2$  test. The stage and location of pressure ulcers were also described and compared. The incidence of pressure ulcers at 31 days was described by the same parameters in the general population described previously.

##### SECONDARY CRITERIA

The stage (NPUAP) of pressure ulcers, their location, caregiver satisfaction with the use of Axtair Automorpho® Plus supports as well as patients' opinions of bed rest comfort, and any adverse events were described by numbers and frequencies. The evolution of the Braden score in the incidence population was described and compared by Student's tests on repeated series.

##### SOFTWARE USED

The data was recorded on Clinsight and analyzed with SAS version 9.3 software in the department of medical biostatistics of Cen Biotech. The significance level was set at  $P < 0.05$ .

##### SAMPLE SIZE

The latest national PERSE survey of 21,538 patients evaluated the current prevalence of pressure ulcers in rehabilitation patients at 11.8%. The number of cases required to compare it in a bilateral situation, with an alpha risk of 5% and a power of 80%, to the expected rate of 7.0% of pressure ulcers when using Axtair Automorpho® plus, according to the "Historical control studies" design is 290 patients<sup>10</sup> rounded up to 330 to allow for 10% of patient drop-out.

##### Information to individuals and informed consent

In line with the legislation, the consent of the person or of a trusted third party was obtained after the practitioner had given them an explanatory letter and answered their questions. People under 18 years of age or with an impaired cognitive state characterized by an MMSE <15 could not participate in the study. The study was approved by the CCTIRS on 15 June 2016 and authorized by the CNIL (DR-2017-025).

**Prevalence study**

All patients present in each department 60 days and 120 days after the start of the study were included in the prevalence population, respectively 231 and 225 patients, making a total of 456 patients. The objective of this observational prevalence study with historical control was to determine whether a significant decrease in the prevalence of pressure ulcers could be achieved by providing Axtair Automorpho® Plus supports for all patients for whom the caregivers considered it necessary.

These 456 patients, of whom 58.8% were female, were 72.3±15.6 years old: 21.9% of patients under 60 years, 40.1% aged 60 to 80 years, and 37.9% over 80 years old. The risk factors for pressure ulcers included malnutrition (11.4%), urinary and/or fecal incontinence (25.9%), hemiplegia or paraplegia (20.7%), bone fracture (23.5%), obesity (BMI> = 30 kg/m<sup>2</sup>) (23.7%), diabetes (21.9%), peripheral neuropathies (3.1%), and psychotropic medication (33.8%).

The beds of 21.5% (N.=98) of these patients were equipped with Axtair Automorpho® Plus supports and the reasons why the others were not equipped were mainly because caregivers considered they did not require it (80.7%). Other reasons related to contraindications (1.4%) and especially to the fact that the patient was already in the department before the start of the study (15.1%). The supports used for the other patients were static mattresses to help prevent bedsores (74.3%), standard mattresses (22.6%), and other motorized air mattresses (3.1%).

**Primary endpoint: total prevalence of pressure ulcers and historical comparison**

On the day of the prevalence study, 18 of the 456 patients presented pressure ulcers which occurred during the period when Axtair Automorpho® Plus supports were made available to the departments. This corresponds to a prevalence of 3.9% [2.4, 6.2]. This figure is significantly lower than those recorded in the PERSE Study<sup>6</sup> chosen for the historical comparison: 11.8% [10.8; 12.8] P<0.0001 (Table I).

TABLE I.—Comparison of the prevalence of patients with pressure ulcers to the prevalence described in the national Perse Study.

Pressure ulcers	CEN 1462		CI 95%	PERSE Study	
	N.	%		N.	%
No	438	96.1	[93.8 ;97.4]	3324	88.2
Yes	18	3.9	[2.4; 6.2]	444	11.8
Total	456	100.0		3768	100.0

χ<sup>2</sup>=25.641, P value: <0.0001; significance: <0.0001.

These 18 patients were 71.2±17.7 years old. Twelve (66.6%) were female and presented a high risk of pressure ulcers with a Braden score on admission to the service of 15.7±3.2. They presented numerous risk factors favoring the occurrence of pressure ulcers and in particular 50.0% presented urinary, fecal or mixed incontinence, 44.4% psychoactive drug intake, 38.9% a paralytic syndrome (hemiplegia or paraplegia), 33.3% a fracture, 27.8% diabetes, 27.8% obesity, and 11.1% under-nutrition. Among them, six patients (33.3%) already had a pressure ulcer.

These 18 patients presented a total of 21 pressure ulcers that occurred despite the prevention policy. These pressure ulcers were located at the heels (52.4%, N.=11), sacrum (28.6%, N.=6), the feet other than the heels (9.5%, N.=2), ischion (4.8%, N.=1), and scrotum (4.8%, N.=1). Of the 11 heel pressure sores, there were 1 stage I, 6 stage II, 2 stage III, and 2 stage IV ulcers. Among the pressure ulcers on the sacrum, there were 2 stage I, 2 stage II, and 2 stage III pressure ulcers. The pressure ulcer on the feet other than the heel were 1 stage II and 1 stage III. The pressure ulcer of the ischium was stage II and that of the scrotum also. A total of 14.3% (N.=3) of pressure ulcers were stage I, 52.4% (N.=11) stage II, 23.8% (N.=5) stage III, and 9.5% (N.=2) stage IV. These stages were not statistically different from those found in the PERSE Study because of the small number of pressure ulcers that reduced the power of the statistical test, but stage II sores were factually more frequent than in the PERSE Study<sup>6</sup> 52.4% vs. 38.3% and conversely stage III and IV sores less frequent 33.3% vs. 44.3%. As for the locations, the same lack of statistical power was present, but it was noticeable that although the locations on the heels occurred with comparable frequency (52.4% vs. 58.6% in the PERSE Study), the locations on the sacrum were less numerous (28.6% vs. 43.3%).

**Prevalence in the moderate and high-risk subpopulation**

An analysis was conducted in moderate to high risk patients targeted, *i.e.* patients without a pressure ulcer on admission to the ward, with a Braden score of 17 or less, bedridden over 15 hours, and up out of bed during the day. In 70 patients in this subpopulation, no pressure ulcers occurred.

**Prevalence in the moderate and very high-risk subpopulation**

In addition, an analysis was carried out in patients with a moderate to very high risk defined by patients with a Braden score of 17 or less, bedridden for more than 15

This document is protected by international copyright laws. No additional reproduction is authorized. It is permitted for personal use to download and save only one file and print only one copy of this Article. It is not permitted to make additional copies (either sporadically or systematically, either printed or electronic) of the Article for any purpose. It is not permitted to distribute the electronic copy of the article through online internet and/or intranet file sharing systems, electronic mailing or any other means which may allow access to the Article. The use of all or any part of the Article for any Commercial Use is not permitted. The creation of derivative works from the Article is not permitted. The production of reprints for personal or commercial use is not permitted. It is not permitted to remove, cover, overlay, obscure, block, or change any copyright notices or terms of use which the Publisher may post on the Article. It is not permitted to frame or use framing techniques to enclose any trademark, logo or other proprietary information of the Publisher.

hours, but who may already have a pressure ulcer and who are not up out of bed during the day. Of the 89 patients concerned, two had a pressure ulcer, *i.e.* a prevalence of 2.2% [0.3; 7.9]. Patient PASO082, was an 83-year-old woman, with hemiplegia, on psychoactive drugs, with a Braden score of 12, and presenting two pressure ulcers; one stage III located on the heel and one stage III on the foot other than the heel. She already had 2 stage II pressure ulcers upon being admitted to the ward, one on the other heel and one not on the heel. Patient PAST015 was a 59 year old man who was quadriplegic, taking psychoactive drugs, and with a Braden score of 11. He had a stage II ulcer located on the heel. He already had one stage II pressure ulcer on the heel of the other foot upon admission to the service.

**Incidence study**

Of the 679 patients admitted to the four wards participating in the study within the five months of inclusion, 57 patients were included in the ulcer incidence study and were hospitalized on a bed with an Axtair Automorpho® Plus support because of their risk of contracting pressure ulcers.

Of these 57 patients, 16 left the study prematurely for the following reasons: four returned home (25%), one was intolerant of the support (6.3%), one presented a contraindication to the support (6.3%), one was considered as no longer requiring prevention with the support (6.3%), seven were transferred to another department (43.8%), and two died (12.5%). Their skin condition was described on the form before leaving the department.

These patients were 79.1±11.7 years old and 57.9% were women. Among their risk factors for pressure ulcers were malnutrition (42.1%), urinary and/or fecal incontinence (38.6%), hemiplegia or paraplegia (33%), bone fracture (26.3%), obesity (BMI≥30 kg/m<sup>2</sup>) (17.5%), diabetes (19.3%), peripheral neuropathy (8.8%), and psychoactive medication (7.0%).

Two-thirds of these patients (66.7%) were bedridden for more than 15 hours a day and one-third were continuously in the chair during the day. Their Braden score was 14.4±3.0 and 82.5% of them had a Braden score of 17 or less.

Of these patients, 40.4% (N.=23) had at least one pressure ulcer and the average number of pressure ulcers presented by the patients with pressure ulcers was 2.5±2.1, ranging from 1 to 9: 1 for 43.5% of patients, 2 for 21.7%, 3 for 13.0%, and 4 and more for 21.7%. In total, these patients had 58 pressure ulcers.

TABLE II.—Description of the pressure sores that occurred during the incidence study.

Patient N.	Heel	Foot other than the heel	Total
PASO126		1 stage II	1
PASO129	1 stage II	2 stage I	3
PAST165		1 stage I	1
PAST2019	1 stage I		1
Total	2	4	

These 58 pressure ulcers were located on the pelvis for 32.8% (N.=19) of them, on the heels for 13.8% (N.=8), on the ischions for 1.7% (N.=1), on the occiput for 1.0% (N.=1). Other localizations included the foot but not the heel for 31.0% (N.=18), the upper limb for 5.2% (N.=3), the lower limb (knee, calf) for 5.2% (N.=3), the genitals (testes and penis) for 3.4% (N.=2), the face for 3.4% (N.=2), and the back for 1.7% (N.=1).

These pressure ulcers were stage I for 27.6% (N.=16), stage II for 43.1% (N.=25), stage III for 15.5% (N.=9), and stage IV for 13.8% (N.=8). It was noted that half of the stage I pressure ulcers were located elsewhere than on the pelvis, heels, and occiput and in particular were located on the foot other than the heel.

One pressure ulcer occurred in 3 patients and 3 pressure ulcers occurred in another patient, *i.e.* 4 out of 57 (Table II) corresponding to an incidence of 7.0% [2,0; 17,0].

These pressure ulcers occurred in:

- patient PASO126, an 83-year-old man, hemiplegic, with mixed incontinence, a Braden score of 13, and a bed rest period >15 hours who presented a stage II ulcer located on the foot other than the heel;
- patient PASO129, an 83-year-old woman, hemiplegic, with mixed incontinence, a Braden score of 12 and bed resting time > 15 hours, who presented 3 ulcers: 2 stage I on the foot other than the heel and 1 stage II ulcer on the heel;
- patient PAST165, a 59-year-old man, quadriplegic, with mixed incontinence, Braden score of 11, and bed rest >15 hours, who had a stage I ulcer on the foot but not the heel;
- patient PAST219, an 80-year-old woman, with Guillain Barré syndrome, with mixed incontinence, Braden score of 13 and bed rest > 15 hours, who presented a stage I pressure ulcer located on the heel.

This patient description highlights that:

- pressure ulcers were stage I for four of them (66.7%) and stage II for the other two (33.3%);
- three patients had stage I pressure ulcers, and two had stage II pressure ulcers, which allows us to evaluate the incidence of stage I pressure ulcers as 3 out of 57 patients, *i.e.* 5.2% [0.1; 14.6], the incidence of stage II pressure ul-

cers as 2 out of 57 patients, *i.e.* 3.5% [0.1, 12.1] and that due to the lack of stage III and IV, the incidence of stages III and IV is 0;

- of the six pressure ulcers that occurred, four occurred on the foot other than the heel (66.0%) and two on the heel, and of the four patients, two presented only localizations on the foot other than the heel and two occurred on the heel, corresponding to an incidence of 3.5% [0.1, 12.1];
- these pressure ulcers all occurred in high-risk patients presenting pathologies with paralysis.

In the subpopulation of reference, namely the 18 patients whose hospital stay exceeded 29 days, who had a Braden score of  $\leq 17$ , and who were bedridden more than 15 hours per day and free of pressure sores upon admission to the department, a pressure ulcer occurred in only 1 of them, *i.e.* an incidence of pressure ulcers of 5.6% [0.1, 27.3]. This pressure ulcer occurred in Patient PASO126, an 83-year-old, male, hemiplegic (right hemiplegia), with mixed incontinence, a Braden score of 13, and bed rest >15 hours and he presented a stage II pressure ulcer located on the right foot other than on the heel.

In addition, the frequency of patients with at least one pressure ulcer fell significantly during the study from 40.4% [27.6, 54.2] at inclusion to 22.4% at D15 [11.8, 36.6] and 19.5% [8.8, 34.9] at D31 (Figure 1).

At day 15, 11 patients had a total of 22 pressure ulcers.


Figure 1.—Patients presenting one or more pressure sores at the different visits.


Figure 2.—Evolution of the Braden Score.  $\Delta J0/J15$ , *t*-test: P value: 0.0003 Significance: <0.001  $\Delta D0/D31$ , *t*-test: P value: 0.0005. Significance: P<0.001.

These 22 pressure ulcers were located on the pelvis for 22.7% (N.=5) of them and 27.3% on the heels (N.=6). The other 11 locations concerned in particular the foot other than the heel for 22.7% (N.=5), the upper limb for 9.1% (N.=2), the lower limb (knee, calf) for 9.1% (N.=2), the face for 4.5% (N.=1), and the back for 4.5% (N.=1). These pressure ulcers were stage I for 31.8% (N.=7) of them, stage II for 40.9% (N.=9), stage III for 13.6% (N.=3), and stage IV for 13.6% (N.=3).

At D31, eight patients still had a total of 17 pressure ulcers located on the heel for 29.4% (N.=5) and on the pelvis for 11.8% (N.=2). As previously, the most frequent locations were the foot other than the heel for 47.1% (N.=8) and the upper limb for 11.8% (N.=2). These pressure ulcers were stage I for 17.6% (N.=3) of them, stage II for 35.3% (N.=6), stage III for 23.5% (N.=4), and stage IV for 23.5% (N.=4).

There was also an improvement in the Braden score throughout the follow-up period, which evolved from 14.4±3.0 at D0 to 15.4±3.1 at D15, and 15.7±3.4 at D31 (<0.001) corresponding to an average increase of 1.3±2.7 points (Figure 2).

The percentages of patients whose skin did not present maceration also increased from 51.1% at inclusion to 63.3% at D15 and 73.2% at D31 (P<0.001) (Figure 3).

The opinions of the patients and of the caregivers were also collected for each patient monitored in the incidence study. More than two-thirds of the patients (70.4%) rated the Axtair Automorpho® Plus mattress comfortable or very comfortable at day 15, and in particular, 90.9% of patients rated the noise level as low or very low.

Axtair Automorpho® Plus supports were also rated very positively by caregivers at D15. All the caregivers found the noise level of the support “low or very low”, 90.9% judged that the support was easy or very easy to set-up and maintain, and 86.4% found it made it very easy to turn patients and sit them up.

In terms of materiovigilance, five patients experienced


Figure 3.—Patients without maceration at the different follow-up visits.

eight adverse events, which were all serious adverse events but none of them were considered to be attributable to the Axtair Automorpho® Plus support. These serious adverse events were five deaths, one of which was preceded by three hospitalizations leading to death. The reasons for death were: subdural hematoma following a fall (patient ETIE037), heart failure (patient HCOL012), acute dehydration (patient HCOL045), and respiratory distress (patients PAST060 and PAST176).

## Results

The results show that no pressure ulcer occurred on the sacrum where the support acts strongly and that two-thirds of the pressure ulcers described in the study were only stage I, which also demonstrates how effective prevention was. They also show that two-thirds of pressure ulcers were located on the foot other than the heel, where the Axtair Automorpho® Plus support is the least likely to act.

This increased awareness also made for more thorough reporting of pressure ulcers during the follow-up of the study. Half of the incidence is represented by patients with stage I sores and localized on the foot other than the heel. This has to be given consideration because most studies attribute little importance to stage I pressure ulcers and ignore pressure ulcers occurring elsewhere than the occiput, sacrum, or heels. If this type of non-exhaustive assessment had been conducted in this study, the incidence of stage II to IV pressure ulcers would have been 3.5%. These elements must be compared with the results of the incidence of stage II to IV pressure ulcers resulting from a systematic review of the literature concerning 13 studies evaluating the effect of prevention with static air mattresses<sup>13</sup> and which shows an average incidence of 7.8% when using such supports. Under these conditions, the ratio of stage II to IV incidences with Axtair Automorpho® Plus and with static support is 0.44 (3.5%/7.8%) and corresponds to the ratio of 0.42 of the incidence of pressure sores with motorized air structures and with classical supports calculated by C. Shi<sup>14</sup> in a meta-analysis of 65 studies.

It is also noticeable that all these pressure sores occurred in severely paralyzed patients. In the subpopulation of high risk patients, the incidence falls to 5.6% [0.1; 27,3] and corresponds to one high risk patient of 83 years old who was hemiplegic, with mixed incontinence, a Braden score of 13, and bed rest of >15 hours per day. In this patient, again, the pressure ulcer was located on the foot other than the heel and no pressure ulcer was described on the occiput, sacrum, or heels.

In addition to information about the rate of pressure ulcer occurrence, the global incidence study also provides information about the evolution of pressure ulcer occurrence in this population of bedridden patients on Axtair Automorpho® Plus. It appears that not only the number of incidents of pressure ulcers is low, but that, far from increasing, the proportion of patients with a pressure ulcer decreases significantly in this population placed on the Axtair Automorpho® Plus support from 40.4% [27.6; 54.2] at inclusion to 22.4% [11.8; 36.6] at day 15, and 19.5% [8.8; 34.9] at day 31, *i.e.* it is halved. During this follow-up period, the percentages of maceration-free patients also increased from 51.1% at inclusion to 63.3% at day 15 and 73.2% at day 31, while the Braden scores rose over the follow-up period. These indications suggest that not only does the patient's use of Axtair Automorpho® Plus prevent the occurrence of pressure ulcers but that it contributes, in conjunction with the other measures implemented in the service and a probable improvement in the patients' health, to an improvement in those sores already present when patients were first put on this support.

## Discussion

The prevalence of 3.9% [2.4, 6.2] observed in the study for all RS patients where Axtair Automorpho® Plus supports were available is statistically lower than that observed for RS services in the PERSE Study<sup>6</sup> chosen for the historical comparison: 11.8% [10.8, 12.8] ( $P < 0.0001$ ) and demonstrates the effectiveness of the prevention afforded by Axtair Automorpho® Plus supports. This quantitative efficiency is reflected qualitatively with 11.0% fewer stage III and stage IV pressure ulcers and 14.7% fewer sacrum locations.

In moderate to very high-risk patients, the prevalence is 2.2% [0.3; 7.9]. This lower prevalence among subjects with a moderate to very high risk than in the general population should not be viewed as paradoxical but as the effect of the prevention tools and of the careful monitoring of the patients at risk. It should be noted that the two patients in whom pressure ulcers occurred were patients at very high risk. One was hemiplegic, the other quadriplegic, taking psychoactive drugs, with Braden scores of 11 or 12 and who already had one or more pressure ulcers and one of them also had mixed incontinence.

The incidence study focuses exclusively on patients whose beds were equipped with Axtair Automorpho® Plus supports because of the risk of pressure ulcers and it shows an incidence of 7.0% [2.0; 17.0]. It is important to underline that in the context of this study, caregivers were

made very much aware of the need to declare stage I ulcers and pressure ulcers occurring at locations other than those most commonly described in literature such as on the foot other than the heel. This is reflected in the results of the study and confirms the evaluation of the preventive effect of the Axtair Automorpho® Plus support.

### Conclusions

The results of the pressure ulcer prevalence study in rehabilitation departments where Axtair Automorpho® Plus supports were available and its comparison with results described in the PERSE national population study in rehabilitation services, as well as results of the study of pressure ulcer incidence in patients bedridden on Axtair Automorpho® Plus, confirm the expected benefit provided to patients in terms of prevention.

### References

1. European Pressure Ulcer Advisory Panel, National Pressure Ulcer Advisory Panel. Prevention and Treatment of Pressure Ulcers: Clinical Practice Guideline, EPUAP, NPUAP, PPPIA. Second edition; 2014 [Internet]. Available from: <https://www.npuap.org/wp-content/uploads/2014/08/Updated-10-16-14-Quick-Reference-Guide-DIGITAL-NPUAP-EPUAP-PPPIA-16Oct2014.pdf> [cited 2018, Jun 7].
2. Conférence de consensus "Prévention et traitement des escarres de

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*Conflicts of interest.*—Brigitte Barrois, Marie-Line Gaubert-Dahan and François-André Allaert received grant to coordinate the study; Ilinca Urzica, Bernard Bordet, Jean-Marc Michel, Marie-Ange Blanchon and Cédric Chol received fees to be investigator of the study.

*Acknowledgments.*—The authors thank the director and paramedical staff of the participating nursing homes in Rehabilitation Services centers: Stéphanie Piot, Stéphanie Grocolas (CRRF COS Pasteur – Troyes); Dominique Delandre, Nikolay Dobrev, Mohammed Yahla, Laurence Leveque, Christelle Tardi, Céline Picard (CRF Pasori – Cosne-sur-Loire); Patricia Blaise, Adeline Rabczynski, Marc Allenbach, Claire Matter, Sophie Mavel (Les Hôpitaux Civils de Colmar – Colmar); Sylvie Favier, Jean-Christophe Poncet (Hôpital Charité du CHU de Saint-Etienne – Saint-Etienne). In addition, the authors thank Winnicare France (manufacturer) and these service providers (1er Mai médical, Bernamont médical, LVM88, HMS42) for providing the medical services needed for this study in accordance with good practice recommendations.

Manuscript accepted: June 4, 2018. - Manuscript received: May 31, 2018.

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KIT COMMERCIAUX	REF SURMATELAS AVEC BASE	REF COMPRESSEUR
VAXT4/AUTO-P VAXT4/CIC-P VAXT4/CIC100-P VAXT4/AUTO100-P VAXT4/AUTO120-P VAXT4/CIC120-P	VAXT/MA/AUTOP VAXT/MA/AUTOP/CIC VAXT/MA/AUTOP100/CIC VAXT/MA/AUTOP100 VAXT/MA/AUTOP120 VAXT/MA/AUTOP120/CIC	VAXT/POMPE/AUTO VAXT/POMPE/AUTO VAXT/POMPE/AUTO VAXT/POMPE/AUTO VAXT/POMPE/AUTO VAXT/POMPE/AUTO

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