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## User manual

### Single mattress for adult medical bed

#### DEVICE INTENDED USE

Compensation for handicap in combination with a medical bed for adult (standard EN 60601-2-52)

#### INDICATIONS

Adult of atypical height - with temporary or permanent loss of autonomy and with a height greater than or equal to 146 cm.

#### CONTRAINDICATIONS

Weight and height of the adult person with atypical morphology, below or above the limits set in the technical data table. For use on medical imaging table and stretcher.

#### ADVERSE SIDE EFFECTS

No adverse effects reported.

**Any serious incident occurring in connection with the device must be notified to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. Inform the competent authority if you consider or have reason to believe that the device poses a serious risk or has been tampered with.**

#### MEDICAL DEVICE COMPOSITION/TECHNICAL DATA

| References Models and versions  | VMB14/HIP  | VMB14/HIPH  |
|---------------------------------|--|---|
| Dimensions (cm)                 | L 200 x l 90 x h 14                                | L 200 x l 90 x 14                                 |
| Foam Specifications             | Polyurethane HR<br>Density<br>34 kg/m <sup>3</sup> | Polyurethane HR<br>Density<br>34kg/m <sup>3</sup> |
| Cover Specification             | Polyurethane / polyester                           | Polyurethane / polyester                          |
| Patient weight min and max (kg) | 30 to 120 kg                                       | 30 to 120 kg                                      |
| Patient height min and max (cm) | 145 <H<185   | 145 <H<185  |

The references mentioned can be completed with a suffix to constitute a commercial reference (example VMB14CC/TH)

#### CLINICAL BENEFIT, DEVICE PERFORMANCE, MECHANISM OF ACTION

DEVICE PERFORMANCE CHARACTERISTICS: This mattress adapts to various positions of the medical bed base necessary to provide treatment, medical procedures, and/or the comfort of bed confinement and basic daily life acts (physical hygiene, transfers, etc.).

Expected clinical benefits: contributes to maintaining help with life activities.

INFORMATION FOR HEALTHCARE PROFESSIONALS: Monitor the condition of the patient's skin in contact with the support. Check the condition of the medical bed base and the protective cover of the mattress foam

#### PREREQUISITES BEFORE USE AND INSTRUCTIONS FOR USE

TRAINING AND QUALIFICATION OF THE DEVICE USER: user training must be carried out by individuals trained and approved by the economic operators concerned, notably in terms of safety and non-compliance reporting.

DEVICE installation: The mattress is delivered with its removable zipped cover. The product is ready to be installed. The cover must be dry.

Preventive maintenance: Carry out a regular visual check of the condition of the foam: visible sagging of material and a non-homogeneous and very slow return of the foam are ageing criteria that compromise the properties of the support. Annual inspection of the condition of the cover (surface aspect and zippers) by exposing the inner surface to a light source: no holes and/or tears. Replace the cover if its surface appearance changes.

#### CLEANING AND DISINFECTION

| Cover |  |  |  |  |  |  |
|-------|---|---|---|---|---|---|
|       | Moderate washing up to 90°C   | Maximum chlorine concentration of 5000ppm   | Do not iron   | Do not dry clean  | Dry with reduced thermal restriction  | Use an approved surface disinfectant/detergent  |

#### INFORMATION ON APPROPRIATE PROCESSES TO ALLOW ITS REUSE

The product must exhibit physical and bacteriological cleanliness. Item treated with a biocide substance at no risk to users.

Avoid direct wet contact with the foam. Prohibit all scouring, stripping products or solvents and "piercing-sharp" objects in direct contact with the cover.

#### WARNINGS, PRECAUTIONS FOR USE, MEASURES REQUIRED

##### PRECAUTIONS FOR USE

In the event of the risk of bedsores, combine with discharge devices or a positioning system or replace the mattress with a therapeutic mattress to help prevent bedsores.

#### WARNING

Assess the absence of risk of the child getting trapped in non mobile parts of the medical bed in accordance with the distances outlined in the standard EN 60601-2-52 and the recommendations set by the manufacturer of the medical bed (read the instruction leaflet). When the user is agitated, the use of full length or segmented lateral barriers with a retractable obstruction accessory must be prescribed. In the event of external cardiac massage, use a board placed between the thorax and the top surface of the mattress.

#### STORAGE, HANDLING, DISPOSAL

##### STORAGE AND USAGE CONDITIONS

The supports must ideally be stored flat, away from direct light and excessive humidity.

|                      | Use                                     | Storage                         |
|----------------------|---|---------------------------------|
| Temperature range    | +15°C to +45°C / +59°F to +113°F        | -25°C to 70°C / +13°F to +158°F |
| Humidity range       | 30% - 70%                               | 30% - 95%                       |
| Atmospheric Pressure | 50kPa - 106 kPa // 7.25 Psi / 15.37 Psi |                                 |

##### SERVICE LIFE

Viscoelastic foam (overlay) mattresses have an estimated service life of 6 years.

##### PRODUCT DISPOSAL

Do not dispose of the product outside of the dedicated area. Respect the established recycling systems in your country.

#### GUARANTEE

The duration of the guarantee for the complete mattress (foam and protective cover) is 2 years. This guarantee starts from the date of purchase of the product from your distributor. This guarantee does not cover normal wear and tear of the product and its protective cover and is not a substitute for legal guarantees. For all complaints, please contact your distributor showing him/her the incriminated product. They will carry out the necessary procedures with our company to proceed with either repair or exchange.

#### Class I medical device according to Regulation (EU) 2017/745.

Technical compliance certificates (FCBA): Passes the fire resistance tests according to the standard EN 597 parts 1 and 2

#### Contents of package

1 mattress  
1 protective cover  
1 user manual

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