WIN CARE

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**DEVICE INTENDED USE**

Viscoelastic foam static mattress, called "memory foam" is intended to be used in combination with a medical bed, for adults or children, for specific medical purposes due to a loss of autonomy and/or deterioration of health status resulting in reduced mobility and the risk of pressure ulcers.

INDICATIONS

For bedridden adults or children patients with mobility impairment with "Medium to High" risk of pressure ulcers (according to the Braden scale or any other validated scale and on medical advice). Assists in the treatment of stage 1 to 2 pressure sores according to EPUAP, up to stage 4 outside the support area.

CONTRAINDICATIONS

Patient weight level below or above the limits defined in the technical data table. For use on medical imaging table.

ADVERSE SIDE EFFECTS

Possible increase in certain sensory impairments in people with multiple sclerosis. Immobility syndrome in the elderly suffering from loss of independence. Instability of seated transfers (excluding mattresses with VMABS stabilizing side borders).

Any serious incident occurring in connection with the device must be notified to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. Inform the competent authority if you consider or have reason to believe that the risk poses a serious risk or has been tampered with.

MEDICAL DEVICE COMPOSITION/TECHNICAL DATA

References Models and versions*	Products for adult medical beds (height ≥ 146 cm)						For Junior children's beds		
	VMA62 VSA65T+ VMS	VMA63 VSA75T+ VMS	VMA64	VMA64S	VMA70 VMA80	VMA117XL ou VMA117XLJRC VMA137XL ou VMA137XLJRC	VMA120BAR VMA140BAR	VMA82-CIC	VMA/PEDIA/M
Sizes (cm)	196x87x14	196x87x14	196x87x14	195x85x14	196x70x14 196x80x14	195x117x17 195x137x17	196x117x17 196x137x17	170x82x14	160x80x14
Specification Foam Sup	PU VE: ≥75kg/m ³ 1.8 kPa	PU VE: ≥75kg/m ³ 1.8 kPa	PU VE: ≥75kg/m ³ 1.8 kPa	PU VE: ≥75kg/m ³ 1.8 kPa	PU VE: ≥75kg/m ³ 1.8 kPa	PU VE: 85 kg/m ³ 1.7 kPa	PU VE: 95kg/m ³ 4.5kPa	PU VE: 85 kg/m ³ 1.7 kPa	PU VE: 85 kg/m ³ 1.7 kPa
Specification Foam Inf	PU: 18kg/m ³ 3.0 kPa	PU: 28kg/m ³ 3.4 kPa	PU: ≥34kg/m ³ 3.4 kPa	PU: ≥34kg/m ³ 3.4 kPa	PU: ≥34kg/m ³ 3.4 kPa	PU: 28 kg/m ³ 3.6 kPa	PU: 40kg/m ³ 5.5 kPa	PU: 28kg/m ³ 3.4 kPa	PU: 28kg/m ³ 3.4 kPa
Patient weight min and max	30 à 120 kg	30 à 150 kg	30 à 180 kg	30 à 180 kg	30 à 135 kg	30 à 150 kg	135 à 270 kg	15 to 70 kg Size: from 125 to 145 cm	15 to 70 kg Size: from 125 to 145 cm
Foam material	100% PU								
Cover material	PES mesh (except CIC with PE mesh) - PU induction								

The references mentioned can be completed with a suffix to constitute a commercial reference according to configuration (type of protective cover, packaging) and/or distributor (examples: VMA62/TH ; VMA64CIC; VMA63BOX-OXP...)

CLINICAL BENEFIT, DEVICE PERFORMANCE, MECHANISM OF ACTION DEVICE PERFORMANCE CHARACTERISTICS: By immersion and enveloping, this support relieves pressure on the body by improving blood circulation and consequently reducing the risk of pressure ulcers in the supporting areas of the body. Integrated heel support.

EXPECTED CLINICAL BENEFITS: Maintenance of tissue oxygenation at the cutaneous and subcutaneous tissue level in contact with the support.

INFORMATION FOR HEALTHCARE PROFESSIONALS: Monitor the condition of the patient's skin in contact with the support several times a day. Use pressure relief devices or positioning systems in patients with pressure ulcer(s).

PREREQUISITES BEFORE USE AND INSTRUCTIONS FOR USE

TRAINING AND QUALIFICATION OF THE DEVICE USER: user training must be carried out by individuals trained and approved by the economic operators concerned, notably in terms of safety and non-compliance reporting.

DEVICE installation: The mattress is delivered with its removable zipped cover. The product is ready to be installed. The surface of the viscoelastic foam must be positioned in contact with the patient's body.

Preventive maintenance: Carry out a regular visual check of the condition of the foam: visible sagging of the ALOVA material and a non-homogeneous and very slow return of the foam are ageing criteria that compromise the properties of the support in helping to prevent pressure sores. Annual inspection of the condition of the cover (surface aspect and zippers) by exposing the inner surface to a light source: no holes and/or tears. Replace the cover if its surface appearance changes.

CLEANING AND DISINFECTION

Cover						
	Moderate washing up to 90°C	Maximum chlorine concentration of 5000ppm	Do not iron	Do not dry clean	Dry with reduced thermal restriction	Use an approved surface disinfectant/detergent

INFORMATION ON APPROPRIATE PROCESSES TO ALLOW ITS REUSE

Before re-using for another patient, the product must exhibit physical and bacteriological cleanliness after application of the above procedures. Do not use the mattress without a cover. Item treated with a biocide substance at no risk to users. www.winncare.fr
Avoid extended wet contact with the foam. Always protect the mattresses with the supplied rigorously dry mattress protectors.
Prohibit all scouring, stripping products or solvents and "piercing-sharp" objects in direct contact with the cover.

WARNINGS, PRECAUTIONS FOR USE, MEASURES REQUIRED**PRECAUTIONS FOR USE**

Non-stabilized bone and/or muscle injuries in contact with the support. COPD-type bronchopulmonary disorders with severe reduction in respiratory muscle tone.

A bedsore is a lesion of varying depth in the skin, associated with excessive and prolonged compression of tissue between the body and a support surface. This excessive pressure can prevent blood circulation and lead to pressure ulcers, a dangerous lesion because it develops from the inside of the body to the outside, and is therefore initially invisible. A pressure sore can have several forms: a simple persistent rash for more than a day, hardening of the skin, a wound of varying depth that can in severe cases affect muscles or underlying bone.

Pressure (max) in mmHg before ageing, at the maximum validated weight, measured flat (0°)

SUITABILITY FOR USE	Max weight User	Head	Buttocks	Heel
VMA Mattress + Promust PU/PUHD/CIC Cover	150 kg	89	60	75
BAR Mattress + Promust PU Cover	270 kg	130	85	93

WARNINGS

(Re) assess the risks of patient entrapment in the non-moving parts of the medical bed associated with "therapeutic mattresses and accessories and articulated positions of the bed base" according to standard EN 60601-2-52 in adults and EN 50637 in adults. child. The use of segmented or full-length side rails should be prescribed.

External cardiac massage not compatible without a board inserted between the thorax and the upper surface of the mattress

STORAGE, HANDLING, DISPOSAL**STORAGE AND USAGE CONDITIONS**

The supports must ideally be stored flat, away from direct light and excessive humidity.

	Use	Storage
Temperature range	+15°C to +45°C / +59°F to +113°F	-25°C to 70°C / +13°F to +158°F
Humidity range	30% - 70%	30% - 95%
Atmospheric pressure	50kPa - 106 kPa // 7.25 Psi / 15.37 Psi	

ROLLED MATTRESSES

ALOVA: 2 mn → 10 mn: This time may be increased if the storage temperature is lower than 15°C and/or in line with the storage period (0 to 6 months). The return to a normal visual appearance of the mattress under normal usage temperature conditions may take up to 24 hours for ALOVA, APLOT and MORFÉA mattresses but this does not present any risk for their use. Small creases may be visible on the lower surface of the CIC cover after unrolling. They present no risk to the patient and do not affect the efficiency of the final product.

SERVICE LIFE

Viscoelastic foam (overlay) mattresses have an estimated service life of 6 years.
The rolled-compressed mattress life span is 6 months (see date on package label).

PRODUCT DISPOSAL

Do not dispose of the product outside of the dedicated area. Respect the established recycling systems in your country.

GUARANTEE

The guarantee duration for ALOVA (overlay) mattresses is 3 years. This guarantee starts from the date of purchase of the product from your distributor. In the case of single patient use, the guarantee relates to the support and its protective cover. In the case of multi-patient use, the duration of guarantee for the protective cover is 2 years. This guarantee does not cover normal wear and tear of the product and its protective cover and is not a substitute for legal guarantees. Please contact your distributor about a defective product. They will carry out the necessary procedures with our company to proceed with either repair or exchange.


Class 1 medical device according to Regulation (EU) 2017/745.

Management systems: Quality: ISO 13485:2016 - Environmental: ISO 14001:2015

Technical compliance certificates (FCBA) : on request or downloadable from our website

Contents of package

- 1 mattress
- 1 protective cover
- 1 user manual

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(max) Pressure in mmHg before ageing, using the maximum validated weight, measured flat (0°)

IDONITA ALLO SCOPO	Max. weight User	Head	Buttocks	Heels
Materasso VMA + Fodera Promust PU/PUHD/CIC	150 kg	89	60	75
Materasso BAR + Fodera Promust PU	270 kg	130	85	93

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The estimated life span of ALOVA mattresses is 6 years.

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GARANZIA

The warranty period for ALOVA mattresses is 3 years. The warranty begins on the purchase date from the distributor. In the case of a single patient use, the warranty covers the support and its protection. In caso di utilizzo su più pazienti, il periodo di garanzia per la protezione è di 2 anni. This warranty does not cover normal wear and tear and is not a substitution of legal guarantees. Please contact your distributor and present the incriminating product. The distributor will carry out the necessary procedure with our company in order to either repair or replace the product.

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