

LOUIS PHILIPPE ROOM



CLASSIC AND REASSURING...




FURNITURE FOR THE LOUIS PHILIPPE ROOM



BOARDS

DIMENSIONS:
Head: W975 x H560mm
Foot: W975 x H700mm


- ▶ Beech veneer and varnished bed boards
- ▶ Widths available: 90, 120, 140 and 160 cm



DESK

DIMENSIONS:
W1100 x D600 x H750 mm


- ▶ 1 drawer
- ▶ Not compatible with Bridge



DRESSER

DIMENSIONS:
W990 x D490 x H850 mm

- ▶ 1 drawer (1)



BEDSIDE CABINET

DIMENSIONS:
W420 x D330 x H640 mm

- ▶ 1 door
- ▶ 1 shelf



BEDSIDE UNIT

DIMENSIONS:
W420 x D320 x H600 mm

- ▶ 1 drawer
- ▶ 1 shelf with no door



TWO-DOOR WARDROBE

DIMENSIONS:
W1360 x D585 x H585 mm

- ▶ Choice of PVC coating
- ▶ 1 side with clothes rail and 1 shelf
- ▶ 1 side with 3 shelves
- ▶ Door opening 110° (invisible hinges)
- ▶ 2 key locks

(1) Drawers are assembled on metal roller runners with built-in end stops, plastic glides with moisture barrier, door opening 110°

LOUIS PHILIPPE	FURNITURE AND BOARD COLORS		SURROUND, EDGING AND HAND RAIL COLORS	MATERIAL
	STANDARD COLORS	OTHER COLORS		
BOARDS	Copper Plum color	-	Copper Plum	Solid beech and re-veneered beech
FURNITURE	Copper Plum color	-		Solid beech and re-veneered beech

1-year warranty*. *Under normal conditions of use and excluding wearing parts - please see instructions of use

Photos not contractually binding

Ambiance Louis Philippe
All in Copper Plum
Legs, molding and banding in solid beech and re-veneered beech
The chair and armchair are no longer available

FP000014/03.2017



MÉDICATLANTIC Groupe WinnCare - 4, Le Pas du Château - 85670 Saint-Paul-Mont-Penit - France
Tel: +33 (0)2 51 98 55 64 - Fax: +33 (0)2 51 98 59 07 - E-mail: info@medicatatlantic.fr - www.winnCare.fr

In the interest of constantly improving its products, the WINNCARE Group reserves the rights to modify the technical specifications of any products without prior notice
Class I medical device in accordance with European Directive 93/42/EEC

