

Assistive devices

- for beach and pool access



WINCARE

N o r d i c

WinnCare Nordic ApS • Taarnborgvej 12C • DK - 4220 Korsoer
Tel. +45 7027 3720 • info@winnCare.dk • www.winnCare.com

Physical and psychological benefits of being in water

Water is a unique environment that provides gentle pressure and soothing support that has many benefits for people with disabilities.

Water exercise has been shown to strengthen muscles in the entire body. Water makes the body lighter and people experience a sense of weightlessness, while gaining a larger range of motion, improved balance, and increased physical endurance.

People with disabilities feel more confident trying movements that they would not attempt on land. This confidence has obvious psychological benefits as well as physical ones.



WinnCare has created two products for easy water access for disabled people; The Marina beach buggy and the Pool E-185 pool lift. If you also need mobility and hygiene solutions for changing rooms and bathrooms as well, do not hesitate to contact us.

Marina beach buggy

Marina is a high quality beach wheelchair that provides wheelchair users an easier way to access and navigate across the beach.

The arm rests are made of floating materials that enable the chair to fully float in water.

The wide wheels can tackle the soft beach sand, making it suitable for all types of terrain.

Marina is made of high strength materials and is quickly disassembled and folded for transport or storage. Customize it with your own logo on the backrest.



youtu.be/B4ke7DQV5pA





Pool E-185

WinnCare's Pool E-185 rotating pool lift allow disabled access to swimming pools, spas and jacuzzis, allowing for a range of hydrotherapy and leisure activities. A solution which is safe and comfortable for the user and easy and effortless for the assistant to use.

Place the user in a sling, a bathing chair or a wide stretcher which all attach easily to the lift.

Product highlights:

- The lift is suitable for people weighing up to 185 kg
- Elevation range up to 2.10 m
- High quality control system from LINAK®
- Emergency stop and powered emergency lowering
- Can be placed from min. 15 cm. to max. 50 cm. from the edge of the pool

See the different transfer equipment and possibilities for mounting below.

Choose the right transfer equipment according to the user's needs and body control:



Bathing chair



Bathing sling




Stretcher

Two solutions for mounting the Pool E-185:




Rotation base to be cast in concrete



Rotation base to be screwed directly into the ground





Do you have questions about one or more products in this brochure or would you like to know more about Winncare Nordic's solutions for patient lifting, hygiene, beds, mattresses and cushions, please feel free to contact us at any time.

Tel.: +45 7027 3720 / E-mail: info@winncare.dk / Web: www.winncare.com

WINCARE
N o r d i c