

BOARDS RANGE

MADÉLIA



- ▶ Light Beech inside finish
- ▶ Ral1001 edging

ABELIA II *



- ▶ Light Beech inside finish
- ▶ Glued ABS edging, Dark Beech color
- ▶ Ral1001 brown epoxy hand rail
- ▶ Ral1001 brown polyethylene ends

AUZENCE II



- ▶ Grey Dragon Tree inside finish
- ▶ Ral 7030 surround

VAL DE VIE *



- ▶ Light Beech inside finish
- ▶ Ral 1001 surround

CARMEN II *



- Headboard**
- ▶ Light Dragon Tree inside finish
 - ▶ matching ABS edging



- Footboard**
- ▶ Light Dragon Tree inside finish
 - ▶ matching ABS edging
 - ▶ Ral 1013 epoxy hand rail

DAGONE *



- ▶ Light Beech inside finish
- ▶ Ral 1001 surround

ELEGANTTI



- ▶ Dune and Chocolate PVC
- ▶ Unique standard colors



CÔTE DU LUMIÈRE *



- ▶ Light Beech inside finish
- ▶ Ral 1001 surround

MEDIDOM II




- Headboard**
- ▶ Light Dragon Tree inside finish
 - ▶ Matching edging



- Footboard**
- ▶ Light Dragon Tree inside finish
 - ▶ Matching edging

PITCHOUNE KALIN



- ▶ Kiwi inside finish
- ▶ Ral 9016 surround

WINNEA®



- Headboard**
- ▶ Crystal Blue finish
 - ▶ Removable headboard



- Footboard**
- ▶ Sixties yellow finish



- Footboard**
- ▶ Crystal Blue finish
 - ▶ Optional blanket holder



Winnea CIC

TECHNICAL SPECIFICATIONS

BOARDS	STANDARD COLORS	BED BASE WIDTH (DIMENSIONS IN CM)						DIMENSIONS: FROM TOP OF BED BASE TO THE TOP OF THE BOARD (HEAD/FOOT)	OVERALL BED DIMENSIONS IN CM (EXCLUDING ACCESSORIES) (DATA TO BE ADDED TO THE BED BASE WIDTH AND LENGTH)		MATERIALS (1)		FIRE RATING
		80	90	100	120	140	160		WIDTH	LENGTH (2)	BOARDS	HAND RAILS	
MADÉLIA II	Light Beech inside finish, Ral 1001 glued ABS edging			x				35	10	6	Chip	-	M3
ABÉLIA II	Light Beech or Noble Cherry inside finish, Ral 1001 glued ABS edging, Ral 8011 brown epoxy hand rail, Ral 8011 brown polyethylene ends		x	x	x	x	x	36	10	7	Chip	Steel	M3
ELEGANTTI	Dune and Chocolate PVC		x					41 / 31	13	10	MDF / PVC / Foam		M1 or M2
CARMEN spaces for wooden side rails			x		x	x	x	46	18	9	MDF / Chip	Steel	M3
CARMEN II	Light Beech inside finish, matching edging and footboard; smooth epoxy hand rail and polyethylene ends in Ral 8011 brown, or Light Dragon Tree inside finish, matching edging and footboard; smooth epoxy hand rail and polyethylene ends in Ral 1013 Cream			x				41	17	9	Chip	Steel (foot)	M3
CARMEN II spaces for wooden side rails	Light Beech inside finish, matching edging and footboard; smooth epoxy hand rail and polyethylene ends in Ral 8011 brown,		x	x	x	x	x	44					
CÔTE DE LUMIÈRE (1)	Light Beech or Noble Cherry inside finish, matching ABS edging uprights and hand rails Ral 1001, lacquer or color		x	x				53	16	10	Solid wood / Chip	Solid wood	M3
VAL DE VIE	Light Beech inside finish, matching ABS edging, Ral 1001 uprights and hand rails		x					46 / 39	6	9	Solid wood / Chip	Solid wood	M3
VAL DE VIE spaces for wooden side rails			x					43,6	13	9			M3
DAGONE	Light Beech or Noble Cherry inside finish, edging and hand rails Dark Beech lacquer		x					49 / 40	11	10	MDF	MDF	M3
AUZENCE II	Grey Dragon Tree interior finish, Ral 7030 taupe lacquer edging Taupe lacquer epoxy hand rails, Ral 7030 upright		x		x	x	x	40	6	11	MDF / Solid wood uprights	Steel	M3
MEDIFORM (1)					x	x	x	53	20	9	MDF / Chip	MDF	M3
MÉDIDOM			x		x	x	x	40 / 13	7	6	MDF	-	M3
MÉDIDOM II	Light Dragon Tree interior finish, matching edging		x	x	x	x	x	37 / 14	6	8	Chip	-	M3
WINNEA FIXES (foot)	Compact Cream Ral 1013 epoxy hand rails			x				42	15	13	Steel / Compact	Steel	M2
WINNEA REMOVABLE (head)								45					
WINNEA CIC	Ral 7047 Grey			x				45	15	13	Polypropylene		M3
PITCHOUNE KALIN	White / Kiwi or White / Violet		x					46	15	6	MDF	-	M3
PITCHOUNE KALIN spaces for wooden side rails			x					43	13	8	MDF	Steel	M3
ACCESSORIES (dimensions in cm)								DIMENSIONS: FROM TOP OF BED BASE TO THE TOP OF THE BOARD					
METAL SIDE RAIL								39	12		-		M2
SAM ACTIV								48	12		-		M3
SAM ERGONOM								46	12		-		M3
SAM EVOLUTION								40	12		-		M3

Our beds and their accessories are designed and manufactured in compliance with their applicable essential requirements as per European Directive 93/42/EEC. We have performed standardized testing on them in their commercial configurations, including the bed boards and accessories that we manufacture, so as to ensure maximum safety and performance for you. As a result, maintenance of the contracted good's warranty depends on the use of original bed boards and accessories, which also guarantees you safe use of the medical bed and its accessories. Médicatlantique shall not be held liable in the event of a combination that is not specified above.

Handset, electric and mechanical parts: 5-year warranty*
Boards and accessories: 1-year warranty*
*Under normal conditions of use and excluding wearing parts - please see instructions of use

(1) boards available with and without spaces for wooden side rails
(2) BED BASE LENGTH: 200 cm for 90/100/120/140/160 cm beds and 160 cm for 80 cm beds
(3) Chip = Chipboard with ABS edging/MDF = medium /Solid = Varnished or lacquered solid wood

Photos not contractually binding

WINCARE Group - 4, Le Pas du Château - 85670 Saint-Paul-Mont-Penit - France


Tél : +33 (0)2 51 98 55 64 - Fax : +33 (0)2 51 98 59 07 - E-mail : info@medicatlantique.fr - www.winnicare.fr

In the interest of constantly improving its products, the WINCARE Group reserves the rights to modify the technical specifications of any product without prior notice

Class I medical device in accordance with European Directive 93/42/EEC



BED ACCESSORIES



DESCRIPTION	ALDRYS' EVOLUTION II							DIMENSIONS	FINISH / MATERIALS
	X'PRESS	X'PRIM 3	AERYS	XXL X'PRESS	XXL DIVISISYS	DUO DIVISISYS			
1 Epoxy side rail	•	•	•	•	•	•	•	Length in top position: 1.47 m Height: 39.5 cm from the bed base	Epoxy
2 Wooden side rail	•	•	•	•	•	•	•	Total bed length Height: 37.5 cm from the bed base	Polymer coating Light Beech
4 Telescopic IV stand 2 plastic hooks	•	•	•	•	•	•	•	Adjustable height: 72 to 141 cm	Stainless steel Plastic hooks
5 Telescopic IV stand 4 plastic hooks	•	•	•	•	•	•	•		
6 IV stand on castors 2 hooks	•	•	•	•	•	•	•	Fixed height: 120 cm	
7 IV stand on castors 4 hooks	•	•	•	•	•	•	•	Adjustable height: 115 to 220 cm	
8 Urinal holder	•	•	•	•	•	•	•		Epoxy
9 Wall-mounted basin holder	•	•	•	•	•	•	•		Epoxy
10 Urine bag hook	•	•	•	•	•	•	•		Epoxy
11 Flexible remote control holder	•	•	•	•	•	•	•		Magnet / Epoxy
12 Wall bumper ALDRYS'EVOLUTION			•						Epoxy
13 Wall bumper X'PRIM 3		•							Epoxy
14 Bed base extension	•	•	•	•	•	•	•	Length: 20 cm	Epoxy
15 Support handle (compatible with side rails and long side panels)	•	•	•	•	•	•	•	Height: 35 cm Length: 29 cm	Epoxy
16 Lifting pole on U base	•	•	•	•	•	•	•	Height: 198 cm	Epoxy
17 Angled lifting pole	•	•	•	•	•	•	•		Epoxy

2-year warranty on accessories*. *Under normal conditions of use and excluding wearing parts - please see instructions of use

Photos not contractually binding