

ÉVILENCE® - CONFORLENCE®



ÉVILENCE table, 1 top
Light Beech color



Conforlence table, 2 tops
Light Dragon Tree color

HYDRAULIC VARIABLE HEIGHT

CONFORLENCE TABLE



2 VERSIONS

- ▶ Single non-tilt top or 2 separate tops.
- ▶ Variable height assisted by gas spring.
- ▶ 19 mm MDF tops, with rounded edges.
- ▶ Top rail on 3 sides to prevent objects from falling off.
- ▶ Top can be tilted into 5 position angles by manual crank.
- ▶ Reading edging strip.

CONVENIENT

- ▶ 11 cm leg height for easy passage under beds.
- ▶ 4 x Ø 50 mm twin-wheel castors, including 2 with brakes.

SAFETY

- ▶ 71.5 cm wide legs for improved stability.

ÉVILENCE TABLE

- ▶ Melamine-faced top, 19 mm thick.
- ▶ Rounded corners.
- ▶ 3-sided top rail
- ▶ Reading edging strip.
- ▶ Table top color: Light Beech.
- ▶ RAL 7035 grey ABS edging (standard color).

2 VERSIONS

- ▶ Single fixed top or 2 separate tops including one tilt top.
- ▶ Option: Two castors with brakes.



FOR HOME CARE

- ▶ The ÉVILENCE table for home use.

FOR USE IN COLLECTIVE CARE

- ▶ The CONFORLENCE table for use in collective care.

ROBUST

- ▶ Steel or aluminum structure, guaranteed for 1 to 2 years.

STYLISH

- ▶ A wide choice of colors to match your interiors.

TECHNICAL SPECIFICATIONS

DESCRIPTION	1 TOP			2 TOPS		RAISING RANGE	LEG HEIGHT	MAX LOAD	CASTORS	WEIGHT
	SINGLE FIXED TOP	SMALL FIXED TOP	LARGE TILT TOP							
ÉVILENCE	92 x 46 cm	-	-	76 to 115 cm	10 cm	15 kg	Ø 50 mm with no brake	13 kg		
	-	29 x 46 cm	64 x 46 cm	76 to 115 cm	10 cm	15 kg	Ø 50 mm with no brake	13 kg		
CONFORLENCE	91 x 43 cm	-	-	76 to 114 cm	11 cm	15 kg	Ø 50 mm with brake	16 kg		
	-	30 x 43 cm	61 x 43 cm	76 to 114 cm	11 cm	15 kg	Ø 50 mm with brake	16 kg		

ÉVILENCE warranty: 1 year - CONFORLENCE warranty: 2 years

Photos not contractually binding

FP0000027/05.2017



WINNCARE Group - 4, Le Pas du Château - 85670 Saint-Paul-Mont-Penit - France
 Tél : +33 (0)2 51 98 55 64 - Fax : +33 (0)2 51 98 59 07 - E-mail : info@medicatlantique.fr - www.winnicare.fr
 In the interest of constantly improving its products, the WINNCARE Group reserves the rights to modify the technical specifications of any product without prior notice
 Class I medical device in accordance with European Directive 93/42/EEC

