



UNIVERSAL BASIC

Perfect sling with additional support for the thighs and body, with a great seating comfort.

USAGE

- For users with lack of stability of torso.
- For transfers from a sitting or lying position.
- For transfers in mobile and stationary hoists.

FEATURES

- The sling is transferring user from a wheelchair/chair or from lying position from a bed or bath.
- The sling provides perfect support around the thighs and the high back, supporting the upper body and giving the user a feeling of comfort and safety.
- The sling ensures a good upright posture that facilitates the positioning of the user.
- Sling in polyester mesh is suitable for bathing.
- The sling has split leg straps.
- Comes with a location pocket for better sling fitting when sitting in a wheelchair, and guide straps on both sides of the sling to make the transfer easier.

SPECIFICATIONS

SWL: 275 kg

Material:

Polyester and polyester mesh.

Sizes: XS-S-M-L-XL-XXL

ITEM NUMBERS:

Universal Basic,
polyester:

	Varenr.
XS	009-09803
S	009-09804
M	009-09805
L	009-09806
XL	009-09807
XXL	009-09808

Universal Basic,
polyester mesh:

	Varenr.
XS	009-09813
S	009-09814
M	009-09815
L	009-09816
XL	009-09817
XXL	009-09818



SLING SIZE GUIDE

The sling has a coloured strap on upper back with size and midline mark, which ensures correct placement of the sling.

Colour code is presented on the size chart.

The image below shows a yellow colour code, which means that the sling is in size M (Medium):



	XS	S	M	L	XL	XXL
Y (hip measurement)	37 cm	42 cm	48 cm	52 cm	68 cm	78 cm
X (measurement from neck to lower back)	50 cm	55 cm	66 cm	75 cm	75 cm	75 cm



Please consult the technical manual for further specifications. All figures and measurements are for guidance only and may differ in actual situations. WinnCare Nordic reserves the right to changes without notice. Product photos may include accessories. Please consult your dealer for further information.



FOR MORE INFORMATION VISIT

winnCare.com

E-mail: info@winnCare.dk

Phone: +45 70273720

WIN  CARE